

# **PATIENT MONITORING SYSTEM USING PIC 16F877**

Copyrights : Mahesh Bharath , 2005

# CONTENTS

## **CONTENTS**

<b>1.</b>	<b>SYNOPSIS</b>	<b>-</b>
<b>2.</b>	<b>BLOCK DIAGRAM</b>	<b>-</b>
<b>3.</b>	<b>INTRODUCTION</b>	<b>-</b>
<b>4.</b>	<b>CIRCUIT DIAGRAM</b>	<b>-</b>
<b>5.</b>	<b>CIRCUIT DESCRIPTION</b>	<b>-</b>
<b>6.</b>	<b>HARDWARE DETAILS</b>	<b>-</b>
<b>7.</b>	<b>SOFTWARE DETAILS</b>	<b>-</b>
<b>8.</b>	<b>CONCLUSION</b>	<b>-</b>
<b>9.</b>	<b>FUTURE ENHANCEMENTS</b>	<b>-</b>
<b>10.</b>	<b>COST ESTIMATE</b>	<b>-</b>
<b>11.</b>	<b>BIBLIOGRAPHY</b>	<b>-</b>

# SYNOPSIS


## **SYNOPSIS**

Our project is a working model which incorporates sensors to measure parameters like body temperature, heart beat rate, respiratory temperature and ECG; and transfer it to the computer so that the patient's health condition can be analyzed by doctors in any part of the hospital. Thus, it reduces the doctor's workload and also gives accurate results. Further this system uses WAP technology which enables the viewing of all parameters on the mobile phone. A micro-controller board is used for analyzing the inputs from the patient and any abnormality felt by the patient causes the monitoring system to give an alarm. Also all the process parameters within an interval selectable by the user are recorded online. This is very useful for future analysis and review of patient's health condition.

For more versatile medical applications, this project can be improvised, by incorporating blood pressure monitoring systems, dental sensors and annunciation systems, thereby making it useful in hospitals as a very efficient and dedicated patient care system.

# **BLOCK DIAGRAM**

## BLOCK DIAGRAM



# INTRODUCTION



## **INTRODUCTION**

### **Measurement of Body Temperature and Respiratory Temperature:**

Thermistor is used for the measurement of body temperature and respiratory temperature. This Thermistor is a passive transducer and its resistance depends on the heat being applied on it. We have arranged the thermistor in the potential divider circuit. This thermistor exhibits a large change in resistance with a change in the body temperature. The respiratory temperature is determined by holding the thermistor near the nose. Initially the thermistor is calibrated to normal body temperature. The thermistor part is attached to the patient whose temperature has to be measured, which changes the resistance value and thus the corresponding change in the temperature is displayed on the monitor graphically. Also all temperature measurements are updated in the patients database. Here in our project we use bead thermistor.

### **Electro Cardiogram (ECG):**

This is a four lead ECG monitoring system. Four sensors are kept at various parts of the body ; the two arms and the two legs. All the signal outputs from the sensors are conditioned by an external circuit consisting of an instrumentation amplifier and a trimpot and is given to the PC through a PIC controller. An interactive program in Visual Basic is developed to read the voltage signals and display a waveform pattern. The four leads used are of silver electrode.

### **Saline Monitoring System:**

For saline monitoring the infrared emitter and detector are placed in a position such that the saline bottle passes between them. They are placed near the neck of the saline bottle.

As long as saline is present, the path of the infrared rays is blocked and the infrared detector is blocked from collecting infrared rays from the infrared emitter. And so the output will indicate normal saline status. The software is written to give an audio alert when the saline level falls below the safe level.

**Patient Calling System:**

The patient calling system consists of four switches which when pressed gives display on the screen and activates an audio alert indicating that a patient is calling. These switches are placed in the vicinity of the patient to enable medical access in an emergency.

**Heart Beat Monitor:**

The patient's heart beat rate is monitored using a photoelectric sensor which can sense the patient's pulse rate. This method of tracking the heart rate is more efficient than the traditional method which derives the same from the ECG graph.


**Attendance Card Reader:**


In order to facilitate easy maintenance and verification of the patient's medical record due knowledge is required regarding the staff interacting with the patient. This is achieved by using a punched identity card and its associated card reader. The card reader uses photo-electric sensors to recognise 16 different medical ward members.



**WAP and Audio Alert:**

Whenever the temperature or the heart beat rate exceed the normal range of values; whenever the saline falls below the preset level and when any patient presses the emergency call switch, provision has been made for an audio alert and also short text messaging the concerned doctor through Wireless Application Protocol.

# **CIRCUIT DIAGRAM**







# **CIRCUIT DESCRIPTION**

**1. Power Supply Circuit:**

All the electronic components only work with a DC supply ranging from  $-12\text{V}$  to  $+12\text{V}$ . We are using the commonly available energy source of  $230\text{V}$ - $50\text{Hz}$  and stepping down, rectifying, filtering and regulating the voltage.

**Step Down Transformer:**

When AC is applied to the primary winding of the power transformer it can either be stepped down or up depending on the value of DC needed. In our circuit the transformer of  $230\text{V}/12-0-12\text{V}$  is used to perform the step down operation where a  $230\text{V}$  AC appears as  $12\text{V}$  AC across the secondary winding. The current rating of the transformer used in our project is  $250\text{mA}$ . Apart from stepping down AC voltages, it gives isolation between the power source and power supply circuitries.

**Rectifier Unit:**

In the power supply unit, rectification is normally achieved using a solid state diode. Diode has the property that will let the electron flow easily in one direction at proper biasing condition. As AC is applied to the diode, electrons only flow when the anode and cathode is negative. Reversing the polarity of voltage will not permit electron flow.

A commonly used circuit for supplying large amounts of DC power is the bridge rectifier. A bridge rectifier of four diodes ( $4 \times \text{IN}4007$ ) are used to achieve full wave rectification. Two diodes will conduct during the negative cycle and the other two will conduct during the positive half cycle. The DC voltage appearing across the output terminals of the bridge rectifier will be somewhat less than 90% of the applied rms value. Normally one



alteration of the input voltage will reverse the polarities. Opposite ends of the transformer will therefore always be 180 deg out of phase with each other.

For a positive cycle, two diodes are connected to the positive voltage at the top winding and only one diode conducts. At the same time one of the other two diodes conducts for the negative voltage that is applied from the bottom winding due to the forward bias for that diode. In this circuit due to positive half cycle D1 & D2 will conduct to give 10.8v pulsating DC. The DC output has a ripple frequency of 100Hz. Since each alteration produces a resulting output pulse, frequency =  $2 \times 50$  Hz. The output obtained is not a pure DC and therefore filtration has to be done.

#### **Filtering Unit:**

Filter circuits which usually capacitor is acting as a surge arrester always follow the rectifier unit. This capacitor is also called as a decoupling capacitor or a bypassing capacitor, is used not only to 'short' the ripple with frequency of 100Hz to ground but also to leave the frequency of the DC to appear at the output. A load resistor R1 is connected so that a reference to the ground is maintained.

1000 $\mu$ f/25v : for the reduction of ripples from the pulsating.

10 $\mu$ f/63v : for maintaining the stability of the voltage at the load side.

0.01 $\mu$ f : for bypassing the high frequency disturbances.

#### **Voltage Regulators:**

The primary purpose of a regulator is to aid the rectifier and filter circuit in providing a constant DC voltage to the device. Power supplies without regulators have an inherent problem of changing DC voltage values due to variations in the load or due to fluctuations in the AC line voltage. With a regulator connected to the DC output, the voltage can be

maintained within a close tolerant region of the desired output. IC7805 and IC7905 is used in this project for providing +5v and -5v DC supply.

## **2. ECG Section:**

The muscular contractions necessary to maintain the hearts' pumping action are initiated by depolarization and repolarization of the SA node and then the depolarization and subsequent repolarization of the AV node. This electrical activity of the heart generates external action potentials. Hence by measuring these potentials at the surface of the body the electrical activity of the heart can be analyzed and the technique used for this is called Electrocardio Graphy.

This is achieved by applying electrodes to certain positions on the body and record the potentials generated between various combinations of these electrodes with an amplifier and the computer screen display. We have used 4 electrodes for the Left Arm, Left Leg, Right Arm and Right Leg. Silver electrodes are preferred as they do not irritate the body and also the contact impedance of Silver is Minimum.

The output from the ECG electrodes are fed to a signal conditioning circuit. This circuit uses high gain amplifiers OP07 to achieve required signed gain. The gain of amplifier is set to meet the minimum input signed requirement of the PLC microcontroller. The OP07 belongs to the family of differential amplifiers which can very efficiently reject the interference of signals which are of a non biological origin. The offset voltage is overcome by using a 500K $\Omega$  trim pot connected for offset balancing of the OP07. The output from each of the four amplifiers are connected to the analog channel inputs AN0-AN3 of the

microcontroller. The controller is programmed to process the ECG signals from a minimum of 3 electrodes.

### 3. Heart Beat Rate Sensor Circuit:

Monitoring the heart beat rate of the patient can be easily accomplished by analyzing the ECG pulse. Here, the ECG pulse is amplified and the average time interval or the instantaneous time interval between two successive R peaks is measured, from which the heart beat rate is derived. But this method fails to indicate heart blocks immediately and so photo electric pulse transducers are used.

The pulse rate monitoring method indicates a heart block immediately by sensing the cessation of blood circulation in the limb terminals. This technique uses photoelectric transducers which are easy to apply than the 3 ECG electrodes. Also the output signal amplitude is large with better signal to noise ratio.

The finger probe used for pulse pick up consists of a Ga As infrared LED and a silicon NPN phototransistor mounted in an enclosure that fits over the tip of the patients' finger. The peak spectral emission of the LED is at 0.94 mm with a 0.707 peak bandwidth of 0.04mm. The silicon phototransistor is sensitive to radiation between 0.4 and 1.1mm. Due to the narrow band of the spectrum involved the radiation heat output is minimized. The photo transistor is used as an emitter follower configuration. The IR signal from the LED is transmitted through the finger tip of the patient's finger and the conductivity of the phototransistor depends on the amount of radiation reaching it with each contraction of the heart, blood is forced to the extremities and amount of blood in finger increases. This alters the optical density and so the IR signal transmission through the finger reduces, causing a correspondence variation in phototransistor output. The phototransistor is connected as part of a voltage divider circuit, with 10K $\Omega$  and 22 K $\Omega$  carbon resistors and produces a voltage

pulse that closely follows the heart beat rate. This pulse output is given to the bit 4 of the Port D of the microcontroller for signal processing.

#### **4. Saline Status Monitoring Circuit:**

The Saline water injection plays a key role in the treatment and recovery of many a patient that it requires constant monitoring. This condition can be easily fulfilled by using IR sensors which can detect a drop in the saline below the quantity. By means of annunciation systems, the hospital staff can be informed and an action of replacing the saline can be easily accomplished before the bottle becomes empty. Also the usage of WAP facilities sending of the saline status to the doctor concerned for any further action required.

The circuit uses an IR emitter and an IR detector which are placed in a straight line with the saline bottle in between, at the point representing the preset saline level. The presence of saline water, in a full bottle, refracts the emitted radiation, thus generating no output at the IR detector. When the saline level falls below the preset value ; the emitted IR radiation causes an photoelectric current output from the detector. The detector output is an analog quantity which is made to drive a switching NPN transistor BC107 to get a binary output from the collector of the transistor. This digital output is fed to the pin 23 of the PIC micro controller, corresponding to Port bit 4 . The signal is processed and the saline status is displayed on the screen. In case of the saline becoming empty the annunciation systems are activated.

#### **5. Attendance Card Reader:**

A patient's medical record is never complete without the data regarding the doctors who treated him/her. In our bid to automate the entire patient care system, this requirement is

also accomplished by using an unique attendance card reader for the medical staff interacting with the patient.

This device is placed along side the patient in his/her ward. It uses an IR emitter- IR detected circuit to read a medical attendant's punched card. The project uses a 4 hole card with four emitter – detector pairs. Thus, a maximum  $2^4 = 16$  staff can be identified using this circuit. Each staff is allotted his punched card and the details of each of the 16 staff is uploaded into the patient's database. Thus, details regarding the date and time of a medical attendant's visit, along with his identity is available for any further reference.

The circuit uses IR detectors which give analog outputs according to the card they read. A set of NPN switching transistors BC107 translate the analog output of the photoelectric card reader in to binary digital signals. These signals are connected to the bits 0 to 3 of Port C (i.e. pin 15 to pin 18) of the PIC micro controller.

#### **6. Patient Call Switches Circuit:**

The patient calling system consists of four switches when pressed gives display on the screen and activates an audio alert indicating that a patient is calling. These switches are placed in the vicinity of the patient to enable medical access in an emergency.

#### **7. Body Temperature Measurement Circuit:**

The temperature measuring circuit uses a thermistor each for the body temperature and the respiratory temperature.

A thermistor is a ceramic semiconductor which exhibits a large change in resistance with a change in its body temperature. The thermistors have much better sensitivity than RTD's and are therefore better suited for precision temperature measurements. The availability of high resistance values allows the thermistors to be used with long extension

leads since the lead resistance or contact resistance effects can be greatly diminished. The non-linearity of the thermistor resistance-temperature characteristics puts a practical limit on the temperature span over which a thermistor can be operated in measurement or control circuit. RTD's have lower sensitivity and are more linear and can therefore be used in applications, where the temperature spans are very wide. Thermistors has other important advantages over RTD's in that they are available in smaller sizes, with faster response times, at lower costs and with greater resistance to shock and vibration effects.

In this circuit we have arranged thermistor in the form of potential divider when thermistor is R1 and a potentiometer is acting as a R2 which forms potential divider network and produces an output from potential divider network which is given to analog input channel of the microcontroller.

In general to obtain clear and constant o/p with respect to the input change, the sensor must be low power consumer. If we draw a lowest current sensitivity the thermistor will improve and provides better performance. Due to the above grounds we have constructed the thermistor circuits to produce low milli volts which can be easily digitized by the PIC. If not the thermistor will to drive large o/p voltage may cause self heating of the device. Self heating means large current flows through the thermistor create heat on it without accepting the body temperature.

#### **8. Respiratory Temperature Measurement Circuit:**

The air is warmed during its passage through the lungs and the expiratory tract and hence there is a detectable difference of temperature between the inhaled air and exhaled air. This difference in the temperature of the exhaled air is the respiratory temperature and it can be sensed by using a thermistor placed in front of the nostrils by means of a suitably holding

device. In case the different of temperature of outside air and that of exhaled air is small, the thermistor can be initially heated to an appropriate temperature and the variation of its resistance in synchronism with the respiration rate, as a result of the cooling effect of the air stream can be detected.

The thermistors with dissipations of about 5 mw to 25 m w are used. It is placed as put of a voltage dividing circuit with a  $330\Omega$  resistance in series with the +5V supply. The output from the voltage divider is given as input to the analog channel of the microcontroller.

#### 9. **PIC Microcontroller Circuit:**

The PIC 16F 877 microcontroller used has a high performance RISC CPU. It is a 40 pin DIP package with many efficient and application friendly features.

We have used a 10 MHz external crystal in crystal oscillator mode during two 22pf capacitors. The crystal output is given to pins 13 and 14 of the DIP.

The microcontroller operates over a voltage range from 2.0v to 5.5v. The +5V input is given to pins 11 and 32 as the positive supply for logic and I10pms, and pin 12, 31 are connected to the ground reference of the +EV regulated power supply board.

The PIC 16F877 has 8 analog channels and out of these analog channels AN0- AN3 are used for input from the ECG signal conditioning circuit, channels AN5 and AN6 are used for the temperature measurement circuits.

The port pins RDO-RD3 are used by the patient call switch circuit and the port pins RCO-RC3 are used by the attendance card reader circuit.

The output from the heart beat sensor is given to bit 4 of port D and output from the saline status monitoring circuit is given to bit 4 of Port C.

Pin number 25 is used in USART asynchronous transmit mode and is connected to the pin 11 of the MAX232 driver.

Pin 26 is used in USART asynchronous receive mode and is connected to pin 12 of MAX 232 IC.

#### **10. MAX 232:**

The MAX 232 power supply section has 2 charge pumps, the first uses external capacitors C1 to double the +5V input to +10V with input impedance of approximately  $200\Omega$ . The second charge pump uses external capacitor to invert +10V to -10V with an overall output impedance of  $45\Omega$ .

The best circuit uses  $22\mu\text{F}$  capacitors for C1 and C4 but the value is not critical. Normally these capacitors are low cost aluminium electrolyte capacitors or tantalum if size is critical. Increasing the value of C1 and C2 to  $47\mu\text{F}$  will lower the output impedance of +5V to +10V doubler by about  $5\Omega$  and +10V to -10V inverter by about  $10\Omega$ . Increasing the value of C3 and C4 lowers the ripple on the power supplies thereby lowering the 16KHz ripple on the RS 232 output. The value of C1 and C4 can be lowered to  $1\mu\text{F}$  in systems where size is critical at the expense of an additional  $20\Omega$  impedance +10V output and  $40\Omega$  additional impedance at -10V input.

#### **Transmitter Section:**

Each of the two transmitters is a CMOS inverter powered by + 10V internally generated supply. The input is TTL and CMOS compatible with a logic threshold of about 26% of  $V_{cc}$ . The input of an unused transmitter section can be left unconnected: an internal  $400\text{K}\Omega$  pull up resistor connected between the transistor input and  $V_{cc}$  will pull the input high forming the unused transistor output low.



The open circuit output voltage swing is guaranteed to meet the RS232 specification +5V output swing under the worst of both transmitter driving the  $3K\Omega$ . Minimum load impedance, the  $V_{cc}$  input at 4.5V and maximum allowable ambient temperature typical voltage with  $5K\Omega$  and  $V_{cc} = +9V$ . The slew rate at output is limited to less than  $30V/\mu s$  and the powered down output impedance will be a minimum of  $300\Omega$  with +2V applied to the output with  $V_{cc} = 0V$ . The outputs are short circuit protected and can be short circuited to ground indefinitely.

### **Receiver Section:**

The two receivers fully conform to RS 232 specifications. Their input impedance is between  $3K\Omega$  either with or without 5V power applied and their switching threshold is within the +3V of RS232 specification. To ensure compatibility with either RS232 or TTL/CMOS input. The MAX232 receivers have  $V_{IL}$  of 0.8V and  $V_{IH}$  of 2.4V the receivers have 0.5V of hysteresis to improve noise rejection. The TTL/CMOS compatible output of receiver will be low whenever the RS232 input is greater than 2.4V. The receiver output will be high when input is floating or driven between +0.8V and -30V.

### **11. RS 232:**

The most common communication interface for short distance is RS-232. RS-232 defines serial communication for one device to one computer communication port, with speeds up to 19,200 baud. Typically 7 or 8 bits (on/off) signal are transmitted to represent a character or digit. 9 pin connector is used.

PIN NAME	25 PIN	9 PIN	I/O DIRECTION
TxD	2	3	Output ("O")
RxD	3	2	Input ("I")
Gnd	7	5	
RTS	4	7	O
CTS	5	8	I
DTR	20	4	O
DSR	6	6	I
RI	22	9	I
DCD	8	1	I

# **HARDWARE DETAILS**



# PIC16F87X

## 28/40-Pin 8-Bit CMOS FLASH Microcontrollers


### Devices Included in this Data Sheet:

- PIC16F873      • PIC16F876
- PIC16F874      • PIC16F877

### Microcontroller Core Features:

- High performance RISC CPU
- Only 35 single word instructions to learn
- All single cycle instructions except for program branches which are two cycle
- Operating speed: DC - 20 MHz clock input  
DC - 200 ns instruction cycle
- Up to 8K x 14 words of FLASH Program Memory,  
Up to 368 x 8 bytes of Data Memory (RAM)  
Up to 256 x 8 bytes of EEPROM Data Memory
- Pinout compatible to the PIC16C73B/74B/76/77
- Interrupt capability (up to 14 sources)
- Eight level deep hardware stack
- Direct, indirect and relative addressing modes
- Power-on Reset (POR)
- Power-up Timer (PWRT) and  
Oscillator Start-up Timer (OST)
- Watchdog Timer (WDT) with its own on-chip RC  
oscillator for reliable operation
- Programmable code protection
- Power saving SLEEP mode
- Selectable oscillator options
- Low power, high speed CMOS FLASH/EEPROM  
technology
- Fully static design
- In-Circuit Serial Programming™ (ICSP) via two  
pins
- Single 5V In-Circuit Serial Programming capability
- In-Circuit Debugging via two pins
- Processor read/write access to program memory
- Wide operating voltage range: 2.0V to 5.5V
- High Sink/Source Current: 25 mA
- Commercial, Industrial and Extended temperature  
ranges
- Low-power consumption:
  - < 0.8 mA typical @ 3V, 4 MHz
  - 20  $\mu$ A typical @ 3V, 32 kHz
  - < 1  $\mu$ A typical standby current

### Pin Diagram




### Peripheral Features:

- Timer0: 8-bit timer/counter with 8-bit prescaler
- Timer1: 16-bit timer/counter with prescaler,  
can be incremented during SLEEP via external  
crystal/clock
- Timer2: 8-bit timer/counter with 8-bit period  
register, prescaler and postscaler
- Two Capture, Compare, PWM modules
  - Capture is 16-bit, max. resolution is 12.5 ns
  - Compare is 16-bit, max. resolution is 200 ns
  - PWM max. resolution is 10-bit
- 10-bit multi-channel Analog-to-Digital converter
- Synchronous Serial Port (SSP) with SPI™ (Master  
mode) and I<sup>2</sup>C™ (Master/Slave)
- Universal Synchronous Asynchronous Receiver  
Transmitter (USART/SCI) with 9-bit address  
detection
- Parallel Slave Port (PSP) 8-bits wide, with  
external  $\overline{RD}$ ,  $\overline{WR}$  and  $\overline{CS}$  controls (40/44-pin only)
- Brown-out detection circuitry for  
Brown-out Reset (BOR)

## PIC16F87X

FIGURE 1-2: PIC16F874 AND PIC16F877 BLOCK DIAGRAM



# PIC16F87X

TABLE 1-2: PIC16F874 AND PIC16F877 PINOUT DESCRIPTION

Pin Name	DIP Pin#	PLCC Pin#	QFP Pin#	I/O/P Type	Buffer Type	Description
OSC1/CLKIN	13	14	30	I	ST/CMOS <sup>(4)</sup>	Oscillator crystal input/external clock source input.
OSC2/CLKOUT	14	15	31	O	—	Oscillator crystal output. Connects to crystal or resonator in crystal oscillator mode. In RC mode, OSC2 pin outputs CLKOUT which has 1/4 the frequency of OSC1, and denotes the instruction cycle rate.
MCLR/VPP	1	2	18	I/P	ST	Master Clear (Reset) input or programming voltage input. This pin is an active low RESET to the device.
RA0/AN0	2	3	19	I/O	TTL	PORTA is a bi-directional I/O port. RA0 can also be analog input0. RA1 can also be analog input1. RA2 can also be analog input2 or negative analog reference voltage. RA3 can also be analog input3 or positive analog reference voltage. RA4 can also be the clock input to the Timer0 timer/counter. Output is open drain type. RA5 can also be analog input4 or the slave select for the synchronous serial port.
RA1/AN1	3	4	20	I/O	TTL	
RA2/AN2/VREF-	4	5	21	I/O	TTL	
RA3/AN3/VREF+	5	6	22	I/O	TTL	
RA4/T0CKI	6	7	23	I/O	ST	
RA5/SS/AN4	7	8	24	I/O	TTL	
RB0/INT	33	36	8	I/O	TTL/ST <sup>(1)</sup>	PORTB is a bi-directional I/O port. PORTB can be software programmed for internal weak pull-up on all inputs. RB0 can also be the external interrupt pin.  RB3 can also be the low voltage programming input. Interrupt-on-change pin. Interrupt-on-change pin. Interrupt-on-change pin or In-Circuit Debugger pin. Serial programming clock. Interrupt-on-change pin or In-Circuit Debugger pin. Serial programming data.
RB1	34	37	9	I/O	TTL	
RB2	35	38	10	I/O	TTL	
RB3/PGM	36	39	11	I/O	TTL	
RB4	37	41	14	I/O	TTL	
RB5	38	42	15	I/O	TTL	
RB6/PGC	39	43	16	I/O	TTL/ST <sup>(2)</sup>	
RB7/PGD	40	44	17	I/O	TTL/ST <sup>(2)</sup>	

Legend: I = input    O = output    I/O = input/output    P = power  
 — = Not used    TTL = TTL input    ST = Schmitt Trigger input

- Note 1: This buffer is a Schmitt Trigger input when configured as an external interrupt.  
 2: This buffer is a Schmitt Trigger input when used in Serial Programming mode.  
 3: This buffer is a Schmitt Trigger input when configured as general purpose I/O and a TTL input when used in the Parallel Slave Port mode (for interfacing to a microprocessor bus).  
 4: This buffer is a Schmitt Trigger input when configured in RC oscillator mode and a CMOS input otherwise.

# PIC16F87X

TABLE 1-2: PIC16F874 AND PIC16F877 PINOUT DESCRIPTION (CONTINUED)

Pin Name	DIP Pin#	PLCC Pin#	QFP Pin#	I/O/P Type	Buffer Type	Description
RC0/T1OSO/T1CKI	15	16	32	I/O	ST	PORTC is a bi-directional I/O port. RC0 can also be the Timer1 oscillator output or a Timer1 clock input. RC1 can also be the Timer1 oscillator input or Capture2 input/Compare2 output/PWM2 output. RC2 can also be the Capture1 input/Compare1 output/PWM1 output. RC3 can also be the synchronous serial clock input/output for both SPI and I <sup>2</sup> C modes. RC4 can also be the SPI Data In (SPI mode) or data I/O (I <sup>2</sup> C mode). RC5 can also be the SPI Data Out (SPI mode). RC6 can also be the USART Asynchronous Transmit or Synchronous Clock. RC7 can also be the USART Asynchronous Receive or Synchronous Data.
RC1/T1OSI/CCP2	16	18	35	I/O	ST	
RC2/CCP1	17	19	36	I/O	ST	
RC3/SCK/SCL	18	20	37	I/O	ST	
RC4/SDI/SDA	23	25	42	I/O	ST	
RC5/SDO	24	26	43	I/O	ST	
RC6/TX/CK	25	27	44	I/O	ST	
RC7/RX/DT	26	29	1	I/O	ST	
RD0/PSP0	19	21	38	I/O	ST/TTL <sup>(3)</sup>	PORTD is a bi-directional I/O port or parallel slave port when interfacing to a microprocessor bus.
RD1/PSP1	20	22	39	I/O	ST/TTL <sup>(3)</sup>	
RD2/PSP2	21	23	40	I/O	ST/TTL <sup>(3)</sup>	
RD3/PSP3	22	24	41	I/O	ST/TTL <sup>(3)</sup>	
RD4/PSP4	27	30	2	I/O	ST/TTL <sup>(3)</sup>	
RD5/PSP5	28	31	3	I/O	ST/TTL <sup>(3)</sup>	
RD6/PSP6	29	32	4	I/O	ST/TTL <sup>(3)</sup>	
RD7/PSP7	30	33	5	I/O	ST/TTL <sup>(3)</sup>	
RE0/ $\overline{\text{RD}}$ /AN5	8	9	25	I/O	ST/TTL <sup>(3)</sup>	PORTE is a bi-directional I/O port. RE0 can also be read control for the parallel slave port, or analog input5. RE1 can also be write control for the parallel slave port, or analog input6. RE2 can also be select control for the parallel slave port, or analog input7.
RE1/ $\overline{\text{WR}}$ /AN6	9	10	26	I/O	ST/TTL <sup>(3)</sup>	
RE2/ $\overline{\text{CS}}$ /AN7	10	11	27	I/O	ST/TTL <sup>(3)</sup>	
Vss	12,31	13,34	6,29	P	—	Ground reference for logic and I/O pins.
VDD	11,32	12,35	7,28	P	—	Positive supply for logic and I/O pins.
NC	—	1,17,28,40	12,13,33,34	—	—	These pins are not internally connected. These pins should be left unconnected.

Legend: I = input    O = output    I/O = input/output    P = power  
 — = Not used    TTL = TTL input    ST = Schmitt Trigger input

- Note 1: This buffer is a Schmitt Trigger input when configured as an external interrupt.  
 2: This buffer is a Schmitt Trigger input when used in Serial Programming mode.  
 3: This buffer is a Schmitt Trigger input when configured as general purpose I/O and a TTL input when used in the Parallel Slave Port mode (for interfacing to a microprocessor bus).  
 4: This buffer is a Schmitt Trigger input when configured in RC oscillator mode and a CMOS input otherwise.

TABLE 3-1: PORTA FUNCTIONS

Name	Bit#	Buffer	Function
RA0/AN0	bit0	TTL	Input/output or analog input.
RA1/AN1	bit1	TTL	Input/output or analog input.
RA2/AN2	bit2	TTL	Input/output or analog input.
RA3/AN3/VREF	bit3	TTL	Input/output or analog input or VREF.
RA4/T0CKI	bit4	ST	Input/output or external clock input for Timer0. Output is open drain type.
RA5/SS/AN4	bit5	TTL	Input/output or slave select input for synchronous serial port or analog input.

Legend: TTL = TTL input, ST = Schmitt Trigger input

TABLE 3-2: SUMMARY OF REGISTERS ASSOCIATED WITH PORTA

Address	Name	Bit 7	Bit 6	Bit 5	Bit 4	Bit 3	Bit 2	Bit 1	Bit 0	Value on: POR, BOR	Value on all other RESETS
05h	PORTA	—	—	RA5	RA4	RA3	RA2	RA1	RA0	--0x 0000	--0u 0000
85h	TRISA	—	—	PORTA Data Direction Register						--11 1111	--11 1111
9Fh	ADCON1	ADFM	—	—	—	PCFG3	PCFG2	PCFG1	PCFG0	--0- 0000	--0- 0000

Legend: x = unknown, u = unchanged, - = unimplemented locations read as '0'.

Shaded cells are not used by PORTA.

TABLE 3-5: PORTC FUNCTIONS

Name	Bit#	Buffer Type	Function
RC0/T1OSO/T1CKI	bit0	ST	Input/output port pin or Timer1 oscillator output/Timer1 clock input.
RC1/T1OSI/CCP2	bit1	ST	Input/output port pin or Timer1 oscillator input or Capture2 input/Compare2 output/PWM2 output.
RC2/CCP1	bit2	ST	Input/output port pin or Capture1 input/Compare1 output/PWM1 output.
RC3/SCK/SCL	bit3	ST	RC3 can also be the synchronous serial clock for both SPI and I <sup>2</sup> C modes.
RC4/SDI/SDA	bit4	ST	RC4 can also be the SPI Data In (SPI mode) or data I/O (I <sup>2</sup> C mode).
RC5/SDO	bit5	ST	Input/output port pin or Synchronous Serial Port data output.
RC6/TX/CK	bit6	ST	Input/output port pin or USART Asynchronous Transmit or Synchronous Clock.
RC7/RX/DT	bit7	ST	Input/output port pin or USART Asynchronous Receive or Synchronous Data.

Legend: ST = Schmitt Trigger input

TABLE 3-6: SUMMARY OF REGISTERS ASSOCIATED WITH PORTC

Address	Name	Bit 7	Bit 6	Bit 5	Bit 4	Bit 3	Bit 2	Bit 1	Bit 0	Value on: POR, BOR	Value on all other RESETS
07h	PORTC	RC7	RC6	RC5	RC4	RC3	RC2	RC1	RC0	xxxx xxxx	xxxx xxxx
87h	TRISC	PORTC Data Direction Register								1111 1111	1111 1111

Legend: x = unknown, u = unchanged



TABLE 3-7: PORTD FUNCTIONS

Name	Bit#	Buffer Type	Function
RD0/PSP0	bit0	ST/TTL <sup>(1)</sup>	Input/output port pin or parallel slave port bit0.
RD1/PSP1	bit1	ST/TTL <sup>(1)</sup>	Input/output port pin or parallel slave port bit1.
RD2/PSP2	bit2	ST/TTL <sup>(1)</sup>	Input/output port pin or parallel slave port bit2.
RD3/PSP3	bit3	ST/TTL <sup>(1)</sup>	Input/output port pin or parallel slave port bit3.
RD4/PSP4	bit4	ST/TTL <sup>(1)</sup>	Input/output port pin or parallel slave port bit4.
RD5/PSP5	bit5	ST/TTL <sup>(1)</sup>	Input/output port pin or parallel slave port bit5.
RD6/PSP6	bit6	ST/TTL <sup>(1)</sup>	Input/output port pin or parallel slave port bit6.
RD7/PSP7	bit7	ST/TTL <sup>(1)</sup>	Input/output port pin or parallel slave port bit7.

Legend: ST = Schmitt Trigger input, TTL = TTL input

Note 1: Input buffers are Schmitt Triggers when in I/O mode and TTL buffers when in Parallel Slave Port mode.

TABLE 3-9: PORTE FUNCTIONS

Name	Bit#	Buffer Type	Function
RE0/ $\overline{RD}$ /AN5	bit0	ST/TTL <sup>(1)</sup>	I/O port pin or read control input in Parallel Slave Port mode or analog input: RD 1 = Idle 0 = Read operation. Contents of PORTD register are output to PORTD I/O pins (if chip selected)
RE1/ $\overline{WR}$ /AN6	bit1	ST/TTL <sup>(1)</sup>	I/O port pin or write control input in Parallel Slave Port mode or analog input: WR 1 = Idle 0 = Write operation. Value of PORTD I/O pins is latched into PORTD register (if chip selected)
RE2/ $\overline{CS}$ /AN7	bit2	ST/TTL <sup>(1)</sup>	I/O port pin or chip select control input in Parallel Slave Port mode or analog input: CS 1 = Device is not selected 0 = Device is selected

Legend: ST = Schmitt Trigger input, TTL = TTL input

Note 1: Input buffers are Schmitt Triggers when in I/O mode and TTL buffers when in Parallel Slave Port mode.

TABLE 3-10: SUMMARY OF REGISTERS ASSOCIATED WITH PORTE

Address	Name	Bit 7	Bit 6	Bit 5	Bit 4	Bit 3	Bit 2	Bit 1	Bit 0	Value on: POR, BOR	Value on all other RESETS
09h	PORTE	—	—	—	—	—	RE2	RE1	RE0	--- -xxx	--- -uuu
89h	TRISE	IBF	OBf	IBOV	PSPMODE	—	PORTE Data Direction Bits			0000 -111	0000 -111
9Fh	ADCON1	ADFM	—	—	—	PCFG3	PCFG2	PCFG1	PCFG0	--0- 0000	--0- 0000

Legend: x = unknown, u = unchanged, - = unimplemented, read as '0'. Shaded cells are not used by PORTE.

TABLE 3-8: SUMMARY OF REGISTERS ASSOCIATED WITH PORTD

Address	Name	Bit 7	Bit 6	Bit 5	Bit 4	Bit 3	Bit 2	Bit 1	Bit 0	Value on: POR, BOR	Value on all other RESETS
08h	PORTD	RD7	RD6	RD5	RD4	RD3	RD2	RD1	RD0	xxxx xxxx	uuuu uuuu
88h	TRISD	PORTD Data Direction Register								1111 1111	1111 1111
89h	TRISE	IBF	OBf	IBOV	PSPMODE	—	PORTE Data Direction Bits			0000 -111	0000 -111

Legend: x = unknown, u = unchanged, - = unimplemented, read as '0'. Shaded cells are not used by PORTD.

**STATUS REGISTER (ADDRESS 03h, 83h, 103h, 183h)**

R/W-0	R/W-0	R/W-0	R-1	R-1	R/W-x	R/W-x	R/W-x
IRP	RP1	RP0	$\overline{TO}$	$\overline{PD}$	Z	DC	C
bit 7							bit 0

bit 7	<b>IRP</b> : Register Bank Select bit (used for indirect addressing) 1 = Bank 2, 3 (100h - 1FFh) 0 = Bank 0, 1 (00h - FFh)	bit 3	<b>PD</b> : Power-down bit 1 = After power-up or by the CLRNDT instruction 0 = By execution of the SLEEP instruction
bit 6-5	<b>RP1:RP0</b> : Register Bank Select bits (used for direct addressing) 11 = Bank 3 (180h - 1FFh) 10 = Bank 2 (100h - 17Fh) 01 = Bank 1 (80h - FFh) 00 = Bank 0 (00h - 7Fh) Each bank is 128 bytes	bit 2	<b>Z</b> : Zero bit 1 = The result of an arithmetic or logic operation is zero 0 = The result of an arithmetic or logic operation is not zero
bit 4	<b>TO</b> : Time-out bit 1 = After power-up, CLRNDT instruction, or SLEEP instruction 0 = A WDT time-out occurred	bit 1	<b>DC</b> : Digit carry/borrow bit (ADDF, ADDLW, SUBLW, SUBWF instructions) (for borrow, the polarity is reversed) 1 = A carry-out from the 4th low order bit of the result occurred 0 = No carry-out from the 4th low order bit of the result
		bit 0	<b>C</b> : Carry/borrow bit (ADDF, ADDLW, SUBLW, SUBWF instructions) 1 = A carry-out from the Most Significant bit of the result occurred 0 = No carry-out from the Most Significant bit of the result occurred

**Note:** For borrow, the polarity is reversed. A subtraction is executed by adding the two's complement of the second operand. For rotate (RRF, RLF) instructions, this bit is loaded with either the high, or low order bit of the source register.

#### INTCON REGISTER (ADDRESS 0Bh, 8Bh, 10Bh, 18Bh)

R/W-0	R/W-0	R/W-0	R/W-0	R/W-0	R/W-0	R/W-0	R/W-x
GIE	PEIE	TOIE	INTE	RBIE	TOIF	INTF	RBIF
bit 7							bit 0


bit 7	<b>GIE</b> : Global Interrupt Enable bit 1 = Enables all unmasked interrupts 0 = Disables all interrupts	bit 3	<b>RBIE</b> : RB Port Change Interrupt Enable bit 1 = Enables the RB port change interrupt 0 = Disables the RB port change interrupt
bit 6	<b>PEIE</b> : Peripheral Interrupt Enable bit 1 = Enables all unmasked peripheral interrupts 0 = Disables all peripheral interrupts	bit 2	<b>TOIF</b> : TMRO Overflow Interrupt Flag bit 1 = TMRO register has overflowed (must be cleared in software) 0 = TMRO register did not overflow
bit 5	<b>TOIE</b> : TMRO Overflow Interrupt Enable bit 1 = Enables the TMRO interrupt 0 = Disables the TMRO interrupt	bit 1	<b>INTF</b> : RBO/INT External Interrupt Flag bit 1 = The RBO/INT external interrupt occurred (must be cleared in software) 0 = The RBO/INT external interrupt did not occur
bit 4	<b>INTE</b> : RBO/INT External Interrupt Enable bit 1 = Enables the RBO/INT external interrupt 0 = Disables the RBO/INT external interrupt	bit 0	<b>RBIF</b> : RB Port Change Interrupt Flag bit 1 = At least one of the RB7:RB4 pins changed state; a mismatch condition will continue the bit. Reading PORTB will end the mismatch condition and allow the bit to be cleared (must be cleared in software). 0 = None of the RB7:RB4 pins have changed state

#### Indirect Addressing, INDF and

##### FSR Registers

The INDF register is not a physical register. Addressing the INDF register will cause indirect addressing. Indirect addressing is possible by using the INDF register. Any instruction using the INDF register actually accesses the register pointed to by the File Select Register, FSR. Reading the INDF register itself, indirectly (FSR = '0') will read 00h. Writing to the INDF register indirectly results in a no operation (although status bits may be affected). An effective 9-bit address is obtained by concatenating the 8-bit FSR register and the IRP bit (STATUS<7>), as shown in figure 2.6.

FIGURE 2-6: DIRECT/INDIRECT ADDRESSING



Note 1: For register file map detail, see Figure 2-3.

## **ANALOG TO DIGITAL CONVERSION**

The Analog-to-Digital (A/D) Converter module has five inputs for the 28-pin devices and eight for the other devices. The analog input charges a sample and hold capacitor. The output of the sample and hold capacitor is the input into the converter. The converter then generates a digital result of this analog level via successive approximation. The A/D conversion of the analog input signal results in a corresponding 10-bit digital number. The A/D module has high and low voltage reference input that is software selectable to some combination of VDD, Vss, RA2, or RA3.

The A/D converter has a unique feature of being able to operate while the device is in SLEEP mode. To operate in SLEEP, the A/D clock must be derived from the

A/D's internal RC oscillator.

The A/D module has four registers. These registers are:

- A/D Result High Register (ADRESH)
- A/D Result Low Register (ADRESL)
- A/D Control Register0 (ADCON0)
- A/D Control Register1 (ADCON1)

The ADCON0 register, shown in Register 11-1, controls the operation of the A/D module. The ADCON1 register, configures the functions of the port pins. The port pins can be configured as analog inputs (RA3 can also be the voltage reference), or as digital I/O. Additional information on using the A/D module can be found in the PICmicro™ Mid-Range MCU Family reference Manual (DS33023)

These steps should be followed for doing an A/D Conversion:

1. Configure the A/D module:
  - Configure analog pins/voltage reference and digital I/O (ADCON1)
  - Select A/D input channel (ADCON0)
  - Select A/D conversion clock (ADCON0)
  - Turn on A/D module (ADCON0)
2. Configure A/D interrupt (if desired):
  - Clear ADIF bit
  - Set ADIE bit
  - Set PEIE bit
  - Set GIE bit
3. Wait the required acquisition time.
4. Start conversion:
  - Set GO/DONE bit (ADCON0)
5. Wait for A/D conversion to complete, by either:
  - Polling for the GO/DONE bit to be cleared (with interrupts enabled); OR
  - Waiting for the A/D interrupt
6. Read A/D result register pair (ADRESH:ADRESL), clear bit ADIF if required.
7. For the next conversion, go to step 1 or step 2, as required. The A/D conversion time per bit is defined as TAD. A minimum wait of 2TAD is required before the next acquisition starts.

### **11.1 A/D Acquisition Requirements**

For the A/D converter to meet its specified accuracy, the charge holding capacitor (CHOLD) must be allowed to fully charge to the input channel voltage level. The analog input model is shown in Figure 11-2. The source impedance (RS) and the internal sampling switch (RSS) impedance directly affect the time required to charge the capacitor CHOLD. The sampling switch (RSS) impedance varies over the device voltage (VDD), see Figure 11-2. The maximum recommended impedance for analog sources is 10 kΩ. As the impedance is decreased, the acquisition time may be decreased.

After the analog input channel is selected (changed), this acquisition must be done before the conversion can be started.

To calculate the minimum acquisition time, Equation 11-1 may be used. This equation assumes that 1/2 LSb error is used (1024 steps for the A/D). The 1/2 LSb error is the maximum error allowed for the A/D to meet its specified resolution.

To calculate the minimum acquisition time, TACQ, see the PICmicro™ Mid-Range Reference Manual (DS33023).

#### **11.4.1 A/D RESULT REGISTERS**

The ADRESH:ADRESL register pair is the location where the 10-bit A/D result is loaded at the completion of the A/D conversion. This register pair is 16-bits wide. The A/D module gives the flexibility to left or right justify

the 10-bit result in the 16-bit result register. The A/D Format Select bit (ADFM) controls this justification. Figure 11-4 shows the operation of the A/D result justification. The extra bits are loaded with '0's'. When an A/D result will not overwrite these locations (A/D disable), these registers may be used as two general purpose 8-bit registers .

### 10.1 USART Baud Rate Generator (BRG)

The BRG supports both the Asynchronous and Synchronous modes of the USART. It is a dedicated 8-bit baud rate generator. The SPBRG register controls the period of a free running 8-bit timer. In Asynchronous mode, bit BRGH (TXSTA<2>) also controls the baud rate. In Synchronous mode, bit BRGH is ignored. Table 10-1 shows the formula for computation of the baud rate for different USART modes which only apply in Master mode (internal clock).

Given the desired baud rate and FOSC, the nearest integer value for the SPBRG register can be calculated using the formula in Table 10-1. From this, the error in baud rate can be determined.

It may be advantageous to use the high baud rate (BRGH = 1), even for slower baud clocks. This is because the  $FOSC/(16(X + 1))$  equation can reduce the baud rate error in some cases.

Writing a new value to the SPBRG register causes the BRG timer to be reset (or cleared). This ensures the BRG does not wait for a timer overflow before outputting the new baud rate.

#### 10.1.1 SAMPLING

The data on the RC7/RX/DT pin is sampled three times by a majority detect circuit to determine if a high or a low level is present at the RX pin.

TABLE 10-1: BAUD RATE FORMULA

SYNC	BRGH = 0 (Low Speed)	BRGH = 1 (High Speed)
0	(Asynchronous) Baud Rate = $FOSC/(64(X+1))$	Baud Rate = $FOSC/(16(X+1))$
1	(Synchronous) Baud Rate = $FOSC/(4(X+1))$	N/A

X = value in SPBRG (0 to 255)

When setting up an Asynchronous Transmission, follow these steps:

1. Initialize the SPBRG register for the appropriate baud rate. If a high speed baud rate is desired, set bit BRGH (Section 10.1).	5. Enable the transmission by setting bit TXEN, which will also set bit TXIF.
2. Enable the asynchronous serial port by clearing bit SYNC and setting bit SPEN.	6. If 9-bit transmission is selected, the ninth bit should be loaded in bit TX9D.
3. If interrupts are desired, then set enable bit TXIE.	7. Load data to the TXREG register (starts transmission).
4. If 9-bit transmission is desired, then set transmit bit TX9.	8. If using interrupts, ensure that GIE and PEIE (bits 7 and 6) of the INTCON register are set.

When setting up an Asynchronous Reception, follow these steps:

1. Initialize the SPBRG register for the appropriate baud rate. If a high speed baud rate is desired, set bit BRGH (Section 10.1).	6. Flag bit RCIF will be set when reception is complete and an interrupt will be generated if enable bit RCIE is set.
2. Enable the asynchronous serial port by clearing bit SYNC and setting bit SPEN.	7. Read the RCSTA register to get the ninth bit (if enabled) and determine if any error occurred during reception.
3. If interrupts are desired, then set enable bit RCIE.	8. Read the 8-bit received data by reading the RCREG register.
4. If 9-bit reception is desired, then set bit RX9.	9. If any error occurred, clear the error by clearing enable bit CREN.
5. Enable the reception by setting bit CREN.	10. If using interrupts, ensure that GIE and PEIE (bits 7 and 6) of the INTCON register are set

TABLE 10-6: REGISTERS ASSOCIATED WITH ASYNCHRONOUS RECEPTION

Address	Name	Bit 7	Bit 6	Bit 5	Bit 4	Bit 3	Bit 2	Bit 1	Bit 0	Value on: POR, BOR	Value on all other RESETS
0Bh, 8Bh, 10Bh, 18Bh	INTCON	GIE	PEIE	TOIE	INTE	RBIE	TOIF	INTF	ROIF	0000 000x	0000 000u
0Ch	PIR1	PSPIF <sup>(1)</sup>	ADIF	RCIF	TXIF	SSPIF	CCP1IF	TMR2IF	TMR1IF	0000 0000	0000 0000
18h	RCSTA	SPEN	RX9	SREN	CREN	—	FERR	OERR	RX9D	0000 -00x	0000 -00x
1Ah	RCREG	USART Receive Register								0000 0000	0000 0000
8Ch	PIE1	PSPIE <sup>(1)</sup>	ADIE	RCIE	TXIE	SSPIE	CCP1IE	TMR2IE	TMR1IE	0000 0000	0000 0000
98h	TXSTA	CSRC	TX9	TXEN	SYNC	—	BRGH	TRMT	TX9D	0000 -010	0000 -010
99h	SPBRG	Baud Rate Generator Register								0000 0000	0000 0000

Legend: x = unknown, - = unimplemented locations read as '0'. Shaded cells are not used for asynchronous reception.


**Note 1:** Bits PSPIE and PSPIF are reserved on PIC16F873/876 devices; always maintain these bits clear.

TABLE 10-5: REGISTERS ASSOCIATED WITH ASYNCHRONOUS TRANSMISSION

Address	Name	Bit 7	Bit 6	Bit 5	Bit 4	Bit 3	Bit 2	Bit 1	Bit 0	Value on: POR, BOR	Value on all other RESETS
0Bh, 8Bh, 10Bh, 18Bh	INTCON	GIE	PEIE	TOIE	INTE	RBIE	TOIF	INTF	ROIF	0000 000x	0000 000u
0Ch	PIR1	PSPIF <sup>(1)</sup>	ADIF	RCIF	TXIF	SSPIF	CCP1IF	TMR2IF	TMR1IF	0000 0000	0000 0000
18h	RCSTA	SPEN	RX9	SREN	CREN	—	FERR	OERR	RX9D	0000 -00x	0000 -00x
19h	TXREG	USART Transmit Register								0000 0000	0000 0000
8Ch	PIE1	PSPIE <sup>(1)</sup>	ADIE	RCIE	TXIE	SSPIE	CCP1IE	TMR2IE	TMR1IE	0000 0000	0000 0000
98h	TXSTA	CSRC	TX9	TXEN	SYNC	—	BRGH	TRMT	TX9D	0000 -010	0000 -010
99h	SPBRG	Baud Rate Generator Register								0000 0000	0000 0000

Legend: x = unknown, - = unimplemented locations read as '0'. Shaded cells are not used for asynchronous transmission.


**Note 1:** Bits PSPIE and PSPIF are reserved on the PIC16F873/876; always maintain these bits clear.




# CODI Semiconductor, Inc.

## LINEAR

### LOGIC SYMBOLS AND CONNECTION DIAGRAMS



### PACKAGE OUTLINES



### LINEAR INTEGRATED CIRCUITS 3 - TERMINAL VOLTAGE REGULATORS

	1 AMP Pos. 7805-7824	1 AMP Neg. 7905-7912	.5A Pos. 7805-7824	.1A Pos. 78L05-78L24
Thermal resistance R <sub>J-θ</sub> R <sub>J-ε</sub>	5°C/W 70°C/W		7°C/W 100°C/W	
Power Dissipation	20W		10W	500mW
Temperature range Operating Ambient Operating Junction Storage	+20 to +80°C +20 to +125°C -55 to +125°C		+20 to +75°C +20 to +125°C -40 to +125°C	
Line Regulation (mV) V <sub>o</sub> /V <sub>i</sub> (%)	2.0		1.0	2.0
Load Regulation (mV) V <sub>o</sub> /I <sub>o</sub> (%)	2.0		1.0	
Logic/Connection diagram	1-CD 1	1-CD 2	1-CD 1	1-CD 3
Package	TO-220		TO-91	

TYPE No.	Nominal OUTPUT Voltage Range (V)	OUTPUT Voltage Range (V)	INPUT Voltage Range (V)	Quiescent Current (mA)	Ripple Rejection (dB) MIN	OUTPUT Noise Voltage (μV)	OUTPUT Voltage Drift (mV/°C)
----------	---	--------------------------------	-------------------------------	------------------------------	---------------------------------	---------------------------------	---------------------------------------

## 1 AMP POSITIVE

7805	5.0	4.8-5.2	7.0-35.0	5.3	62	60	-0.4
7806	6.0	5.75-6.25	8.0-35.0	6.7	59	65	-0.4
7808	8.0	7.7-8.3	10.5-35.0	8.7	56	50	-0.4
7812	12.0	11.5-12.5	14.5-35.0	13.7	55	30	-0.8
7815	15.0	14.4-15.6	17.5-35.0	16.7	54	45	-1.0
7818	18.0	17.3-18.7	21.0-35.0	19.7	53	85	-1.3
7824	24.0	23.0-25.0	27.0-60.0	24.0	50	110	-1.6

## 1 AMP NEGATIVE

7905	-5.0	-4.8-5.2	-7.0-35.0	5.3	56	100	-0.6
7906	-6.0	-5.75-6.25	-8.0-35.0	6.7	54	150	-0.5
7908	-8.0	-7.7-8.3	-10.5-35.0	8.7	58	200	-0.6
7912	-12.0	-11.5-12.5	-14.5-35.0	13.7	56	100	-0.8
7915	-15.0	-14.4-15.6	-17.5-35.0	16.7	54	125	-1.0
7918	-18.0	-17.3-18.7	-21.0-35.0	19.7	54	450	-1.0
7924	-24.0	-23.0-25.0	-27.0-60.0	24.0	51	600	-1.0

## .5A (500ma) POSITIVE

7805	5.0	4.8-5.2	7.0-35.0	6.5	62	60	-1.0
7806	6.0	5.75-6.25	8.0-35.0	6.6	59	65	-1.0
7808	8.0	7.7-8.3	10.5-35.0	6.6	56	50	-1.0
7812	12.0	11.5-12.5	14.5-35.0	8.8	55	75	-1.0
7815	15.0	14.4-15.6	17.5-35.0	8.8	54	90	-1.0
7818	18.0	17.3-18.7	21.0-35.0	8.8	53	100	-1.0
7824	24.0	23.0-25.0	27.0-60.0	9.0	50	170	-1.0

## .1A (100ma) POSITIVE

78L05	5.0	4.75-5.25	7.0-30.0	6.0*	60	40	-1.0
78L06	6.0	5.95-6.05	8.0-30.0	6.0*	59	50	-1.0
78L08	8.0	7.6-8.4	10.5-30.0	6.0*	59	60	-1.0
78L12	12.0	11.6-12.4	14.5-30.0	6.5*	57	80	-1.3
78L15	15.0	14.3-15.7	17.5-30.0	6.5*	56	90	-1.3
78L18	18.0	17.1-18.9	21.0-30.0	6.5*	57	120	-2.0
78L24	24.0	22.8-25.2	27.0-60.0	7.0*	52	200	-2.0


\*max





## BT THERMISTOR

The BT thermistor is a small thermal sensing device providing high reliability, stable characteristics and a wide operating range of -50°C to 300°C. It is used in various applications including medical apparatus, industrial equipment and home electric appliances.

Part number



Dimensions



Specifications

Port No.	R <sub>25</sub> <sup>1</sup>	B value <sup>2</sup>	Dissipation factor (mW/°C)	Thermal Time constant <sup>3</sup>	Rated power at 25°C (mW)	Operating temp. range (°C)
1BT-5	1.000kΩ±10%	3,250K±3%	0.5	4~12	2.5	-50~150
2BT-5	2.000kΩ±10%	3,420K±3%	0.5	4~12	2.5	-50~300
5BT-5(6)	5.000kΩ±10%	3,450K±3%	0.5(0.4)	4~12(3~8)	2.5(2)	-50~300
9BT-5(6)	9.000kΩ±10%	3,470K±3%	0.5(0.4)	4~12(3~8)	2.5(2)	-50~150
10BT-5(6)	10.00kΩ±10%	3,250K±3%	0.5(0.4)	4~12(3~8)	2.5(2)	-50~300
20BT-5(6)	20.00kΩ±10%	3,330K±3%	0.5(0.4)	4~12(3~8)	2.5(2)	-50~300
30BT-5(6)	30.00kΩ±10%	3,450K±3%	0.5(0.4)	4~12(3~8)	2.5(2)	-50~300
40BT-5(6)	40.00kΩ±10%	3,550K±3%	0.5(0.4)	4~12(3~8)	2.5(2)	-50~300
100BT-5(6)	100.0kΩ±10%	3,750K±3%	0.5(0.4)	4~12(3~8)	2.5(2)	-50~300
400BT-5(6)	400.0kΩ±10%	4,050K±3%	0.5(0.4)	4~12(3~8)	2.5(2)	-50~300
500BT-5(6)	500.0kΩ±10%	3,760K±3%	0.5(0.4)	4~12(3~8)	2.5(2)	-50~300
1.3MBT-5(6)	1300kΩ±10%	4,380K±3%	0.5(0.4)	4~12(3~8)	2.5(2)	-50~300

<sup>1</sup> R<sub>25</sub>: Rated zero-power resistance value at 25°C, ±5% are also available.


<sup>2</sup> B value: determined by rated zero-power resistance at 25°C and 85°C.

<sup>3</sup> Time when thermistor temperature reaches 63.2% of the temperature difference. The value is measured in the air.

Resistance-Temperature

Temperature (°C)	Type												
	1BT	2BT	5BT	9BT	10BT	20BT	30BT	40BT	100BT	400BT	500BT	1.3MBT	
-50	29.87	70.71		344.8	284.8	686.7			4860				
-40	17.33	40.02		190.3	163.4	383.1			2599	11043			
-30	10.35	23.28		109.1	97.62	222.6			1439	6198			
-20	6.374	13.96	35.44	64.81	60.41	134.2		306.8	827.4	3573			
-10	4.038	8.640	21.84	39.70	38.63	83.61		184.7	491.1	2109		8066	
0	2.629	5.513	13.87	25.15	25.45	53.75	83.80	115.0	301.4	1274	1486	4598	
10	1.755	3.610	9.057	16.36	17.22	35.53	54.46	73.88	190.1	788.4	945.3	2718	
20	1.200	2.421	6.060	10.92	11.92	24.09	36.37	48.77	123.1	498.7	614.7	1652	
25	1.000	2.000	5.000	9.000	10.00	20.00	30.00	40.00	100.0	400.0	500.0	1300	
30	0.8380	1.661	4.148	7.458	8.434	16.70	24.88	33.00	81.71	322.4	408.9	1029	
40	0.5973	1.163	2.898	5.200	6.084	11.81	17.39	22.82	55.39	212.8	277.7	656.0	
50	0.4338	0.8311	2.065	3.698	4.458	8.511	12.40	16.10	38.31	143.3	192.1	427.8	
60	0.3205	0.6043	1.497	2.677	3.303	6.248	8.990	11.57	27.00	98.24	135.3	284.5	
70	0.2407	0.4468	1.104	1.970	2.460	4.658	6.629	8.457	19.38	68.52	98.90	193.0	
80	0.1834	0.3357	0.8267	1.473	1.850	3.522	4.962	6.279	14.14	46.50	70.47	133.0	
90	0.1417	0.2559	0.6280	1.117	1.405	2.698	3.767	4.730	10.48	34.91	52.00	92.76	
100	0.1110	0.1978	0.4836	0.8581	1.078	2.093	2.897	3.611	7.866	25.47	38.90	65.53	
110		0.08789	0.1547	0.3771	0.6685	0.8355	1.640	2.255	5.968	18.83	29.49	46.91	
120			0.1224	0.2975	0.5264	0.6540	1.296	1.774	4.580	14.10	22.63	34.03	
130			0.09789	0.2373	0.4191	0.5171	1.027	1.410	3.551	10.68	17.55	25.02	
140			0.07908	0.1913	0.3369	0.4126	0.8190	1.130	2.782	8.180	13.76	18.59	
150			0.06450	0.1556	0.2735	0.3321	0.6581	0.9113	1.107	2.201	6.332	10.90	
160				0.1278		0.2696	0.5327	0.7374	0.8943	1.757	4.957	8.717	
170				0.1058		0.2207	0.4345	0.5897	0.7252	1.416	3.915	7.034	
180				0.08836		0.1818	0.3567	0.4909	0.5915	1.150	3.120	5.722	
190				0.07432		0.1508	0.2947	0.4048	0.4859	0.9418	2.508	4.692	
200				0.06295		0.1258	0.2451	0.3360	0.4018	0.7770	2.032	3.876	
210						0.1056	0.2052	0.2806	0.3344	0.6458	1.658	3.225	
220								0.2360	0.2802	0.5403	1.363	2.702	
230								0.1995	0.2381	0.4551	1.127	2.277	
240								0.1695	0.2000	0.3655	0.9390	1.930	
250								0.1447	0.1704	0.3286	0.7864	1.644	
260								0.1243	0.1458	0.2816	0.6626	1.407	
270								0.1072	0.1255	0.2426	0.5613	1.209	
280								0.09283	0.1084	0.2100	0.4780	1.041	
290								0.08078	0.09408	0.1826	0.4091	0.8995	
300								0.07060	0.08200	0.1595	0.3518	0.7810	

Unit(Ω)



February 1995

## LM78XX Series Voltage Regulators

### General Description

The LM78XX series of three terminal regulators is available with several fixed output voltages making them useful in a wide range of applications. One of these is local on card regulation, eliminating the distribution problems associated with single point regulation. The voltages available allow these regulators to be used in logic systems, instrumentation, HiFi, and other solid state electronic equipment. Although designed primarily as fixed voltage regulators these devices can be used with external components to obtain adjustable voltages and currents.

The LM78XX series is available in an aluminum TO-3 package which will allow over 1.0A load current if adequate heat sinking is provided. Current limiting is included to limit the peak output current to a safe value. Safe area protection for the output transistor is provided to limit internal power dissipation. If internal power dissipation becomes too high for the heat sinking provided, the thermal shutdown circuit takes over preventing the IC from overheating.

Considerable effort was expended to make the LM78XX series of regulators easy to use and minimize the number

of external components. It is not necessary to bypass the output, although this does improve transient response. Input bypassing is needed only if the regulator is located far from the filter capacitor of the power supply.

For output voltage other than 5V, 12V and 15V the LM117 series provides an output voltage range from 1.2V to 57V.


### Features

- Output current in excess of 1A
- Internal thermal overload protection
- No external components required
- Output transistor safe area protection
- Internal short circuit current limit
- Available in the aluminum TO-3 package


### Voltage Range

LM7805C	5V
LM7812C	12V
LM7815C	15V

### Schematic and Connection Diagrams




Metal Can Package  
TO-3 (K)  
Aluminum



Bottom View

Order Number LM7805CK,  
LM7812CK or LM7815CK  
See NS Package Number KC02A

Plastic Package  
TO-220 (T)



Top View

Order Number LM7805CT,  
LM7812CT or LM7815CT  
See NS Package Number T03B

**Absolute Maximum Ratings**

If Military/Aerospace specified devices are required, please contact the National Semiconductor Sales Office/Distributors for availability and specifications.

Input Voltage ( $V_O = 5V, 12V$  and  $15V$ ) 35V  
 Internal Power Dissipation (Note 1) Internally Limited  
 Operating Temperature Range ( $T_A$ )  $0^\circ\text{C}$  to  $+70^\circ\text{C}$

Maximum Junction Temperature

(K Package) 150°C  
 (T Package) 150°C  
 Storage Temperature Range  $-65^\circ\text{C}$  to  $+150^\circ\text{C}$   
 Lead Temperature (Soldering, 10 sec.)  
 TO-3 Package K 300°C  
 TO-220 Package T 230°C

**Electrical Characteristics LM78XXC** (Note 2)  $0^\circ\text{C} \leq T_J \leq 125^\circ\text{C}$  unless otherwise noted.

Output Voltage			5V			12V			15V			Unit
Input Voltage (unless otherwise noted)			10V			19V			23V			
Symbol	Parameter	Conditions	Min	Typ	Max	Min	Typ	Max	Min	Typ	Max	
V <sub>O</sub>	Output Voltage	T <sub>J</sub> = 25°C, 5 mA ≤ I <sub>O</sub> ≤ 1 A	4.8	5	5.2	11.5	12	12.5	14.4	15	15.6	V
		P <sub>D</sub> ≤ 15W, 5 mA ≤ I <sub>O</sub> ≤ 1 A	4.75	5.25		11.4	12.6		14.25	15.75		V
		V <sub>MIN</sub> ≤ V <sub>IN</sub> ≤ V <sub>MAX</sub>	(7.5 ≤ V <sub>IN</sub> ≤ 20)			(14.5 ≤ V <sub>IN</sub> ≤ 27)			(17.5 ≤ V <sub>IN</sub> ≤ 30)			V
ΔV <sub>O</sub>	Line Regulation	I <sub>O</sub> = 500 mA, T <sub>J</sub> = 25°C		3	50		4	120		4	150	mV
		ΔV <sub>IN</sub>		(7 ≤ V <sub>IN</sub> ≤ 25)		14.5 ≤ V <sub>IN</sub> ≤ 30)		(17.5 ≤ V <sub>IN</sub> ≤ 30)			V	
		0°C ≤ T <sub>J</sub> ≤ +125°C			50		120		150			mV
		ΔV <sub>IN</sub>		(8 ≤ V <sub>IN</sub> ≤ 20)		(15 ≤ V <sub>IN</sub> ≤ 27)		(18.5 ≤ V <sub>IN</sub> ≤ 30)			V	
		I <sub>O</sub> ≤ 1 A			50		120		150			mV
		ΔV <sub>IN</sub>		(7.5 ≤ V <sub>IN</sub> ≤ 20)		(14.6 ≤ V <sub>IN</sub> ≤ 27)		(17.7 ≤ V <sub>IN</sub> ≤ 30)			V	
ΔV <sub>O</sub>	Load Regulation	T <sub>J</sub> = 25°C			10		12	120		12	150	mV
		5 mA ≤ I <sub>O</sub> ≤ 1.5 A			50		60		75			mV
		250 mA ≤ I <sub>O</sub> ≤ 750 mA			25		60		75			mV
		5 mA ≤ I <sub>O</sub> ≤ 1 A, 0°C ≤ T <sub>J</sub> ≤ +125°C			50		120		150			mV
		I <sub>O</sub>			8		8		8			mA
		0°C ≤ T <sub>J</sub> ≤ +125°C			8.5		8.5		8.5			mA
ΔI <sub>O</sub>	Quiescent Current Change	I <sub>O</sub> ≤ 1 A			0.5		0.5		0.5			mA
		T <sub>J</sub> = 25°C, I <sub>O</sub> ≤ 1 A			1.0		1.0		1.0			mA
		V <sub>MIN</sub> ≤ V <sub>IN</sub> ≤ V <sub>MAX</sub>			(7.5 ≤ V <sub>IN</sub> ≤ 20)		(14.8 ≤ V <sub>IN</sub> ≤ 27)		(17.9 ≤ V <sub>IN</sub> ≤ 30)			V
		I <sub>O</sub> ≤ 500 mA, 0°C ≤ T <sub>J</sub> ≤ +125°C			1.0		1.0		1.0			mA
		V <sub>MIN</sub> ≤ V <sub>IN</sub> ≤ V <sub>MAX</sub>			(7 ≤ V <sub>IN</sub> ≤ 25)		(14.5 ≤ V <sub>IN</sub> ≤ 30)		(17.5 ≤ V <sub>IN</sub> ≤ 30)			V
		V <sub>N</sub>	Output Noise Voltage	T <sub>A</sub> = 25°C, 10 Hz ≤ f ≤ 100 kHz			40		75		90	
$\frac{\Delta V_{IN}}{\Delta V_{OUT}}$	Ripple Rejection	f = 120 Hz { I <sub>O</sub> ≤ 1 A, T <sub>J</sub> = 25°C or I <sub>O</sub> ≤ 500 mA, 0°C ≤ T <sub>J</sub> ≤ +125°C	62	80		55	72		54	70		dB
		V <sub>MIN</sub> ≤ V <sub>IN</sub> ≤ V <sub>MAX</sub>			(8 ≤ V <sub>IN</sub> ≤ 18)		(15 ≤ V <sub>IN</sub> ≤ 25)		(18.5 ≤ V <sub>IN</sub> ≤ 28.5)			V
		I <sub>O</sub> ≤ 1 A, T <sub>J</sub> = 25°C			2.0		2.0		2.0			V
		f = 1 kHz			8		18		19			mΩ
R <sub>O</sub>	Dropout Voltage	T <sub>J</sub> = 25°C, I <sub>OUT</sub> = 1 A			2.0		2.0		2.0			V
		Output Resistance			8		18		19			mΩ
		Short-Circuit Current			2.1		1.5		1.2			A
		Peak Output Current			2.4		2.4		2.4			A
		Average TC of V <sub>OUT</sub>			0.6		1.5		1.8			mV/°C
V <sub>IN</sub>	Input Voltage Required to Maintain Line Regulation	T <sub>J</sub> = 25°C, I <sub>O</sub> ≤ 1 A			7.5		14.6		17.7			V

Note 1: Thermal resistance of the TO-3 package (K, KG) is typically  $4^\circ\text{C/W}$  junction to case and  $35^\circ\text{C/W}$  case to ambient. Thermal resistance of the TO-220 package (T) is typically  $4^\circ\text{C/W}$  junction to case and  $50^\circ\text{C/W}$  case to ambient.

Note 2: All characteristics are measured with capacitor across the input of  $0.22\text{ }\mu\text{F}$ , and a capacitor across the output of  $0.1\text{ }\mu\text{F}$ . All characteristics except no/iso voltage and ripple rejection ratio are measured using pulse techniques ( $t_w \leq 10\text{ ms}$ , duty cycle  $\leq 5\%$ ). Output voltage changes due to changes in internal temperature must be taken into account separately.

MAX232, MAX232I  
DUAL EIA-232 DRIVER/RECEIVER

SLLS047G – FEBRUARY 1989 – REVISED AUGUST 1998


- Operates With Single 5-V Power Supply
- LinBiCMOS™ Process Technology
- Two Drivers and Two Receivers
- ±30-V Input Levels
- Low Supply Current . . . 8 mA Typical
- Meets or Exceeds TIA/EIA-232-F and ITU Recommendation V.28
- Designed to be Interchangeable With Maxim MAX232
- Applications
  - TIA/EIA-232-F
  - Battery-Powered Systems
  - Terminals
  - Modems
  - Computers
- ESD Protection Exceeds 2000 V Per MIL-STD-883, Method 3015
- Package Options Include Plastic Small-Outline (D, DW) Packages and Standard Plastic (N) DIPs

description


The MAX232 device is a dual driver/receiver that includes a capacitive voltage generator to supply EIA-232 voltage levels from a single 5-V supply. Each receiver converts EIA-232 inputs to 5-V TTL/CMOS levels. These receivers have a typical threshold of 1.3 V and a typical hysteresis of 0.5 V, and can accept ±30-V inputs. Each driver converts TTL/CMOS input levels into EIA-232 levels. The driver, receiver, and voltage-generator functions are available as cells in the Texas Instruments LinASIC™ library.

The MAX232 is characterized for operation from 0°C to 70°C. The MAX232I is characterized for operation from –40°C to 85°C.

D, DW, OR N PACKAGE  
(TOP VIEW)



logic symbol†



† This symbol is in accordance with ANSI/IEEE Std 91-1984 and IEC Publication 617-12.

AVAILABLE OPTIONS

T <sub>A</sub>	PACKAGED DEVICES		
	SMALL OUTLINE (D)	SMALL OUTLINE (DW)	PLASTIC DIP (N)
0°C to 70°C	MAX232D‡	MAX232DW‡	MAX232N
–40°C to 85°C	MAX232ID‡	MAX232IDW‡	MAX232IN

‡ This device is available taped and reeled by adding an R to the part number (i.e., MAX232DR).



Please be aware that an important notice concerning availability, standard warranty, and use in critical applications of Texas Instruments semiconductor products and disclaimers thereto appears at the end of this data sheet.

LinASIC and LinBiCMOS are trademarks of Texas Instruments Incorporated.

PRODUCTION DATA Information is current as of publication date. Products conform to specifications per the terms of Texas Instruments standard warranty. Production processing does not necessarily include testing of all parameters.



POST OFFICE BOX 655303 • DALLAS, TEXAS 75265

Copyright © 1998, Texas Instruments Incorporated

## MAX232, MAX232I DUAL EIA-232 DRIVER/RECEIVER

SLLS047G – FEBRUARY 1989 – REVISED AUGUST 1998

### absolute maximum ratings over operating free-air temperature range (unless otherwise noted)<sup>†</sup>

Input supply voltage range, $V_{CC}$ (see Note 1)	–0.3 V to 6 V
Positive output supply voltage range, $V_{S+}$	$V_{CC} - 0.3$ V to 15 V
Negative output supply voltage range, $V_{S-}$	–0.3 V to –15 V
Input voltage range, $V_I$ : Driver	–0.3 V to $V_{CC} + 0.3$ V
Receiver	±30 V
Output voltage range, $V_O$ : T1OUT, T2OUT	$V_{S-} - 0.3$ V to $V_{S+} + 0.3$ V
R1OUT, R2OUT	–0.3 V to $V_{CC} + 0.3$ V
Short-circuit duration: T1OUT, T2OUT	Unlimited
Package thermal impedance, $\theta_{JA}$ (see Note 2): D package	113°C/W
DW package	105°C/W
N package	78°C/W
Storage temperature range, $T_{stg}$	–65°C to 150°C
Lead temperature 1,6 mm (1/16 inch) from case for 10 seconds	260°C

<sup>†</sup> Stresses beyond those listed under "absolute maximum ratings" may cause permanent damage to the device. These are stress ratings only, and functional operation of the device at these or any other conditions beyond those indicated under "recommended operating conditions" is not implied. Exposure to absolute-maximum-rated conditions for extended periods may affect device reliability.

NOTE 1: All voltage values are with respect to network ground terminal.

2. The package thermal impedance is calculated in accordance with JESD 51, except for through-hole packages, which use a trace length of zero.

### recommended operating conditions

	MIN	NOM	MAX	UNIT
Supply voltage, $V_{CC}$	4.5	5	5.5	V
High-level input voltage, $V_{IH}$ (T1IN, T2IN)	2			V
Low-level input voltage, $V_{IL}$ (T1IN, T2IN)			0.8	V
Receiver input voltage, R1IN, R2IN			±30	V
Operating free-air temperature, $T_A$	MAX232	0	70	°C
	MAX232I	–40	85	

# **SOFTWARE DETAILS**

**ASSEMBLY LANGUAGE PROGRAM - PIC16F877**

```

list      p=16f877          ; list directive to define processor
#include <p16f877.inc>       ; processor specific variable definitions

Rrdy_Flag    EQU    0x20
Sync_Flag    EQU    0x21

Count        EQU    0x22
Temp         EQU    0x23

H_Byte       EQU    0x24
L_Byte       EQU    0x25
R0           EQU    0x26      ; RAM Assignments
R1           EQU    0x27
R2           EQU    0x28

H_Temp       EQU    0x29      ; temporary register
L_Temp       EQU    0x2A      ; temporary register

Byte0        EQU    0x2b
Byte1        EQU    0x2c
Byte2        EQU    0x2d
Byte3        EQU    0x2e

Rbuf0        EQU    0x31      ; 0x21 TO 0x28
Rbuf1        EQU    0x31
Rbuf2        EQU    0x32
Rbuf3        EQU    0x33
Rbuf4        EQU    0x34
Rbuf5        EQU    0x35
Rbuf6        EQU    0x36
Rbuf7        EQU    0x37
Rbuf8        EQU    0x38
Rbuf9        EQU    0x39

;*****
; Reset and Interrupt Vectors

    org      00000h ; Reset Vector
    goto     Start

    org      00004h ; Interrupt vector
;*****
; Interrupt Service Routine

IntVector
    ; save context (WREG and STATUS registers) if needed.
    banksel   PIR1
    btfss     PIR1,RCIF    ; Did USART cause interrupt?
    goto      ISREnd

    bcf       PIR1,RCIF
    moviwi    06h          ; mask unwanted bits

```

```

        andwf RCSTA,w      ; check for error
        btfss STATUS,z     ; was error status but set?
        Goto RcvError      ; found error flag it

Receive_Int

        btfss Sync_Flag,0
        goto Syn_Rx
        movlw ']'
        Subwf RCREG,W
        btfss STATUS,Z
        goto Rxchar
        bsf Rrdy_Flag,0
        bcf Sync_Flag,0
        goto ISREnd

Rxchar
        btfss RCREG,6
        goto Numerals
        movlw 0X0F
        andwf RCREG,W
        addlw 0X09

Rxpro
        movwf INDF
        incf FSR,F
        goto ISREnd

Numerals
        movlw 0X0F
        andwf RCREG,W
        goto Rxpro

Syn_Rx
        movlw '['
        subwf RCREG,W
        btfss STATUS,Z
        goto ISREnd
        bsf Sync_Flag,0
        goto ISREnd      ; go to end of ISR, restore context, return

RcvError

        bcf RCSTA,CREN     ; Clear receiver status
        bsf RCSTA,CREN

ISREnd
        ; Restore context if needed.
        retfie

;*****
; Program begins here

        org 00200h ; Beginning of program EPROM
Start
        banksel TRISA

        movlw 0xff
        movwf TRISA

        movlw 0xff
        movwf TRISB

```



```

    movlw 0xbf
    movwf TRISC

    movlw 0xff
    clrf  TRISD

    movlw 0x40    ; 9600 baud @10MHz
    banksel  SPBRG
    movwf  SPBRG

    banksel  PIR1
    bcf  PIR1,RCIF    ; Clear RCIF Interrupt Flag

    banksel  PIE1
    bsf  PIE1,RCIE    ; Set RCIE Interrupt Enable

    banksel  ADCON0
    movlw 0x81
    movwf  ADCON0

    banksel  ADCON1
    movlw 0x80
    movwf  ADCON1

    banksel  INTCON
    bsf  INTCON,PEIE    ; Enable peripheral interrupts
    bsf  INTCON,GIE    ; Enable global interrupts

    bcf  STATUS,RP0
    bcf  STATUS,RP1

    movlw 0x20
    movwf  FSR
    clrf  INDF

ClearNext
    incf  FSR,F
    clrf  INDF
    movf  FSR,W
    xorlw 0x39
    btfss STATUS,Z
    goto  ClearNext

;*****
; Main loop

Wait
    btfss Rrdy_Flag,0 ;s
    goto  Wait

    clrf  Rrdy_Flag

    movf  Rbuf0,W
    xorlw 0X01    ;
    btfss STATUS,Z    ; s
    goto  ChkNext1

    swapf Rbuf1,F    ;
    movf  Rbuf2,W
    iorwf Rbuf1,W

```

```

    andiw 0x10
    ioriw 0xef
    bsf STATUS,RP0 ;
    movwf TRISA ;SET DATA DIRECTION FOR PORT A
    bcf STATUS,RP0 ;

    swapf Rbuf3,F ;
    movf Rbuf4,W ;
    iorwf Rbuf3,W ;WREG = RXBUF3 || RXBUF4
    bsf STATUS,RP0 ;
    movwf TRISB ;SET DATA DIRECTION FOR PORT B
    bcf STATUS,RP0 ;

    swapf Rbuf5,F ;
    mov Rbuf6,w
    iorwf Rbuf5,W ;WREG = RXBUF5 || RXBUF6
    andiw 0xbf
    ioriw 0x80
    bsf STATUS,RP0 ;
    movwf TRISC ;SET DATA DIRECTION FOR PORT C
    bcf STATUS,RP0 ;

    swapf Rbuf7,F ;
    movf Rbuf8,W ;
    iorwf Rbuf7,W ;WREG = RXBUF7 || RXBUF8
    bsf STATUS,RP0 ;
    movwf TRISD ;SET DATA DIRECTION FOR PORT A
    bcf STATUS,RP0 ;SELECT BANK 0

    goto Wait ;

ChkNext1

    movf Rbuf0,W
    xorlw 0x02 ;
    btfss STATUS,Z ; s
    goto ChkNext2

Wt0

    btfss PIR1,TXIF ;s
    goto Wt0
    movlw '{ '
    movwf TXREG

    call Delay
    call Delay

Wt1

    btfss PIR1,TXIF ;s
    goto Wt1

    swapf PORTA,W
    andlw 0x0f
    bsf PCLATH,1
    bsf PCLATH,2
    call HexTable
    bcf PCLATH,1
    bcf PCLATH,2
    movwf TXREG

    call Delay
    call Delay

```

Wt2

```

btfss PIR1,TXIF;s
goto  Wt2

movf   PortA,w

andlw  0x0f
bsf    PCLATH,1
bsf    PCLATH,2
call   HexTable
bcf    PCLATH,1
bcf    PCLATH,2
movwf  TXREG

call   Delay
call   Delay

```

Wt3

```

btfss PIR1,TXIF ;s
goto  Wt3

swapf  PORTB,W
andlw  0x0f
bsf    PCLATH,1
bsf    PCLATH,2
call   HexTable
bcf    PCLATH,1
bcf    PCLATH,2
movwf  TXREG

call   Delay
call   Delay

```

Wt4

```

btfss PIR1,TXIF;s
goto  Wt4

movf   PORTB,W
andlw  0x0f
bsf    PCLATH,1
bsf    PCLATH,2
call   HexTable
bcf    PCLATH,1
bcf    PCLATH,2
movwf  TXREG

call   Delay
call   Delay

```

Wt5

```

btfss PIR1,TXIF ;s
goto  Wt5

swapf  PORTC,W
andlw  0X0F
bsf    PCLATH,1
bsf    PCLATH,2
call   HexTable
bcf    PCLATH,1
bcf    PCLATH,2

```

```

        movwf TXREG

        call Delay
        call Delay

Wt6
        btfss PIR1,TXIF;s
        goto Wt6

        movf PORTC,W
        andlw 0x0f
        bsf PCLATH,1
        bsf PCLATH,2
        call HexTable
        bcf PCLATH,1
        bcf PCLATH,2
        movwf TXREG

        call Delay
        call Delay

Wt7
        btfss PIR1,TXIF ;s

        goto Wt7

        swapf PORTD,W
        andlw 0x0f
        bsf PCLATH,1
        bsf PCLATH,2
        call HexTable
        bcf PCLATH,1
        bcf PCLATH,2
        movwf TXREG

        call Delay
        call Delay

Wt8
        btfss PIR1,TXIF;s
        goto Wt8

        movf PORTD,W
        andlw 0x0f
        bsf PCLATH,1
        bsf PCLATH,2
        call HexTable
        bcf PCLATH,1
        bcf PCLATH,2
        movwf TXREG

        call Delay
        call Delay

Wt9
        btfss PIR1,TXIF
        goto wt9
        movlw '}'
        movwf TXREG

```

```

    call    Delay
    call    Delay

    goto    Wait

ChkNext2

    movf    Rbuf0,W
    xorlw   0X03          ;
    btfss   STATUS,Z      ; s
    goto    ChkNext3

    swapf   Rbuf1,F          ;
    movf    Rbuf2,W
    iorwf   Rbuf1,W
    movwf   PORTA

    swapf   Rbuf3,F          ;
    movf    Rbuf4,W          ;
    iorwf   Rbuf3,W          ;WREG = RXBUF3 || RXBUF4
    movwf   PORTB

    swapf   Rbuf5,F          ;
    movf    Rbuf6,W          ;
    iorwf   Rbuf5,W          ;WREG = RXBUF5 || RXBUF6
    movwf   PORTC

    swapf   Rbuf7,F          ;
    movf    Rbuf8,W          ;
    iorwf   Rbuf7,W          ;WREG = RXBUF7 || RXBUF8
    movwf   PORTD

    goto    Wait          ;

ChkNext3

    movf    Rbuf0,W
    xorlw   0X04          ;
    btfss   STATUS,Z      ; s
    goto    Wait

    rlf     Rbuf1,F
    rlf     Rbuf1,F
    rlf     Rbuf1,F

    movf    Rbuf1,W
    andlw   B'00111000'
    iorlw   0x81

    movwf   ADCON0

    call    Delay          ; Acq Time Delay

    bsf     ADCON0,GO      ;START CONVERSION

WfEoc    btfsc   ADCON0,GO      ;DONE
        goto    WfEoc

    movf    ADRESH,W      ;
    movwf   H_Byte        ;

    bsf     STATUS,RP0
    movf    ADRESL,W      ;

```

```

        bcf     STATUS,RP0
        movwf   L_Byte           ; H_BYTE:L_BYTE = ADRESH:ADRESL

AWt0
        btfss   PIR1,TXIF ;s
        goto    AWt0
        movlw   '{ '
        movwf   TXREG

        call    Delay
        call    Delay

AWt1
        btfss   PIR1,TXIF ;s
        goto    AWt1

        swapf   H_Byte,W
        andlw   0x0f
        bsf     PCLATH,1
        bsf     PCLATH,2
        call    HexTable
        bcf     PCLATH,1
        bcf     PCLATH,2
        movwf   TXREG

        call    Delay
        call    Delay

AWt2
        btfss   PIR1,TXIF;s
        goto    AWt2

        movf    H_Byte,W
        andlw   0x0f
        bsf     PCLATH,1
        bsf     PCLATH,2
        call    HexTable
        bcf     PCLATH,1
        bcf     PCLATH,2
        movwf   TXREG

        call    Delay
        call    Delay

AWt3
        btfss   PIR1,TXIF ;s
        goto    AWt3

        swapf   L_Byte,W
        andlw   0x0f
        bsf     PCLATH,1
        bsf     PCLATH,2
        call    HexTable
        bcf     PCLATH,1
        bcf     PCLATH,2
        movwf   TXREG

        call    Delay
        call    Delay

AWt4
        btfss   PIR1,TXIF;s

```

```

        goto    AWt4

        movf    L_Byte,W
        andlw   0x0f
        bsf     PCLATH,1
        bsf     PCLATH,2
        call    HexTable
        bcf     PCLATH,1
        bcf     PCLATH,2
        movwf   TXREG

        call    Delay
        call    Delay

AWt5:
        btfss   PIR1,TXIF
        goto    AWt5
        movlw   '}'
        movwf   TXREG

        call    Delay
        call    Delay

        goto    Wait

Delay:
        movlw   0xff

movwf   Count

Dwait:
        decfsz  Count,F
        goto    Dwait

        return


        org     0x600

HexTable
        addwf   PCL,f
        retlw   0x30
        retlw   0x31
        retlw   0x32
        retlw   0x33
        retlw   0x34
        retlw   0x35
        retlw   0x36
        retlw   0x37
        retlw   0x38
        retlw   0x39
        retlw   0x41

        retlw   0x42
        retlw   0x43
        retlw   0x44
        retlw   0x45
        retlw   0x46

        end

```

**USER END SOFTWARE IN VISUAL BASIC****Flow Chart:**



**PROGRAM:****Form4.frm**

```

Private Sub BTempTimer_Timer()
MSComm1.Output = "[45]"
Sleep (20)
Rtext = MSComm1.Input
Rtext = Mid$(Rtext, 2, 4)
Pt2 = T2
T2 = Val("&h" & Rtext)
T2 = T2 / 5
T2 = -T2+195
T2=Round(T2,1)

If T2 < 77 Then T2 = Pt2
End If

Text7.text =T2
Label4.Caption = "BodyTemp Vs Time"
DrawWidth = 2
Line (536, 184)-(536, 384), vbYellow
Line (536, 384)-(744, 384), vbYellow
Line (528, 192)-(536, 184), vbYellow
Line (544, 192)-(536, 184), vbYellow
Line (736, 376)-(744, 384), vbYellow
Line (736, 392)-(744, 384), vbYellow
DrawWidth = 2
Line (oldx, oldy)-(536 + k, 384 - (T2 / (200 / 200#))), vbGreen
DrawWidth = 1
oldx = 536 + k
oldy = 384 - (T2 / (200 / 200#))
k = k + 3

If k > 208 Then
Line (528, 182)-(744, 384), Form4.BackColor, BF
oldx = 536#
oldy = 384#
k = 0
End If

End Sub

'Private Sub btimer_Timer()
'btimer.Enabled = False
'out = out And &HCF

```

```
'End Sub
```

```
Private Sub Command1_Click()
MSComm1.Output = "[3ffffbfff]"
Sleep (50)
End
End Sub
```

```
Private Sub Command3_Click()
Form4.Hide
Form3.Show
End Sub
```

```
Private Sub Ecg1Timer_Timer()
MSComm1.Output = "[40]"
Sleep (50)
Rtext = MSComm1.Input
Rtext = Mid$(Rtext, 2, 4)
Ecgl = Val("&h" & Rtext)
Ecgl = Round(Ecgl, 1)
Text9.Text = Ecgl
End Sub
```

```
Private Sub Form_Load()
For i = 1 To 100
Egra1(i) = 50
Next
ptr = 19
MSComm1.PortOpen = True
MSComm1.Output = "[1ffffbfff]"
Sleep (30)
MSComm1.Output = "[3ffffbfff]"
Sarr(1) = 0
Sarr(2) = 0
Sarr(3) = 0
Sarr(4) = 0
Sarr(5) = 0
Sarr(6) = 0
Fnames(1) = "D:\PanimalarCBSS\hitemp.wav"
Fnames(2) = "D:\PanimalarCBSS\Patient1.wav" 'p1 call
Fnames(3) = "D:\PanimalarCBSS\Patient2.wav" 'p2 call
Fnames(4) = "D:\PanimalarCBSS\Patient3.wav" 'p3 call
Fnames(5) = "D:\PanimalarCBSS\Patient4.wav" 'p4 call
Fnames(6) = "D:\PanimalarCBSS\Saline1.wav" 'p4 call
Fcnt = 0
Label5.Caption = Format(Date, "dd/mm/yyyy")
```

```

Set db = OpenDatabase("D:\PanimalarCBSS\hos12.mdb")
Set rs = db.OpenRecordset("patient")
Set Srs = db.OpenRecordset("Digital")
Set Drs = db.OpenRecordset("Doc")
Pa = &H0
Sal = &H0
oldx = 536#
oldy = 384#
k = 0
x=0
a(0) = 5 a(1) = 5 a(2) = 5 a(3) = 6 a(4) = 5 a(5) = 5 a(6) = 5 a(7) = 5 a(8) = 3
a(9) = 10 a(10) = 1 a(11) = 5 a(12) = 5 a(13) = 5 a(14) = 7 a(15) = 5 a(16) = 6
a(17) = 5 a(18) = 5 a(19) = 5

Lg = Shell("C:\Program Files\Microsoft Office\Office\outlook.exe", vbNormalFocus)
End Sub

Private Sub Heartbeattimer_Timer()
MSComm1.Output = "[2]"
Sleep (20)
Rtext = MSComm1.Input
Rtext = Mid$(Rtext, 8, 1)
Beat = Val("&h" & Rtext)
If Beat = 0 And Bflg = 0 Then
Bcnt = Bcnt + 1
Bflg = 1
ElseIf Beat > 0 And Bflg = 1 Then
Bflg = 0
End If
Text1.Text = Round(Bcnt, 0)
End Sub

Private Sub Hrtcnd_Click()
hcnt = 0
Bcnt = 0
End Sub

Private Sub MMControl1_Done(NotifyCode As Integer)
MMControl1.Command = "Close"
Timer2.Enabled = True
End Sub

Private Sub Option1_Click()
Rtemptimer.Enabled = False
Btemptimer.Enabled = False
Ecg1Timer.Enabled = False
Heartbeattimer.Enabled = False

```

```
Rtemptimer.Enabled = True
Picture1.Cls
Text1.Text = Clear
Text5.Text = Clear
Text7.Text = Clear
Text9.Text = Clear
grval = 0
End Sub
```

```
Private Sub Option2_Click()
Rtemptimer.Enabled = False
BTemptimer.Enabled = False
EcglTimer.Enabled = False
Heartbeattimer.Enabled = False
BTemptimer.Enabled = True
Picture1.Cls
Text1.Text = Clear
Text5.Text = Clear
Text7.Text = Clear
Text9.Text = Clear
grval = 1
End Sub
```

```
Private Sub Option3_Click()
Rtemptimer.Enabled = False
BTemptimer.Enabled = False
EcglTimer.Enabled = False
Heartbeattimer.Enabled = False
EcglTimer.Enabled = True
Picture1.Cls
Text1.Text = Clear
Text5.Text = Clear
Text7.Text = Clear
Text9.Text = Clear
End Sub
```

```
Private Sub Option7_Click()
Rtemptimer.Enabled = False
BTemptimer.Enabled = False
EcglTimer.Enabled = False
Heartbeattimer.Enabled = False
Heartbeattimer.Enabled = True
Picture1.Cls
Text1.Text = Clear
Text5.Text = Clear
Text7.Text = Clear
Text9.Text = Clear
```

```

grval = 2
End Sub
Private Sub Rtemptimer_Timer()
MSComm1.Output = "[46]"
Sleep (20)
Rtext = MSComm1.Input
Rtext = Mid$(Rtext, 2, 4)
Pt1 = T1
T1 = Val("&h" & Rtext)
T1 = T1 / 5
T1 = -T1+195
T1 = Round(T1, 1)

If T1 < 77 Then T1 = Pt1
EndIf

Text5.Text = T1
Label4.Caption = "RespTemp Vs Time"
DrawWidth = 2
Line (536, 184)-(536, 384), vbYellow
Line (536, 384)-(744, 384), vbYellow
Line (528, 192)-(536, 184), vbYellow
Line (544, 192)-(536, 184), vbYellow
Line (736, 376)-(745, 385), vbYellow
Line (736, 392)-(744, 384), vbYellow
DrawWidth = 2
Line (oldx, oldy)-(536 + k, 384 - (T1 / (200 / 200#))), vbGreen
DrawWidth = 1
oldx = 536 + k
oldy = 384 - (T1 / (200 / 200#))
k = k + 3
If k > 208 Then
Line (528, 182)-(744, 384), Form4.BackColor, BF
oldx = 536#
oldy = 384#
k = 0
End If

End Sub

Private Sub Timer1_Timer()
Label6.Caption = Format(Time, "hh:mm:ss")

If Heartbeattimer.Enabled = True Then
Text2.Visible = True
Hrtcnd.Visible = True
Else

```

```

Text2.Visible = False
Hrtcnd.Visible = False
End If

If (Ecgl > 100) And Strt = 0 Then
Strt = 1
End If

MSComm1.Output = "[2]"
Sleep (20)
Rtext = MSComm1.Input
Pa = Val("&h" & Mid$(Rtext, 8, 2))
If ((Pa And &H1) = &H0) Then
Label23.Visible = True
p1 = "Yes"
Sarr(2) = 1

If pflg1 = 0 Then
pflg1 = 1
Set objOutlookMsg = objOutlook.CreateItem(olMailItem)
With objOutlookMsg
.To = "9884169070@sms.sancharnet.in" ' The To Address Field
.Subject = "P1 Call" ' txtSubject - The Subject Field
.Body = "Patient 1 is calling"
.Importance = olImportanceHigh ' Priority. Mostly Not Required
.Send ' Send the Composed Message...
End With
Set objOutlookMsg = Nothing
Else
Label23.Visible = False
p1 = "No"
Sarr(2) = 0
End If

If ((Pa And &H2) = &H0) Then
Label24.Visible = True
p2 = "Yes"
Sarr(3) = 1
Else
Label24.Visible = False
p2 = "No"
End If

If ((Pa And &H4) = &H0) Then
Label25.Visible = True
p3 = "Yes"
Sarr(4) = 1

```

```
Else
Label25.Visible = False
p3 = "No"
End If
```

```
If ((Pa And &H8) = &H0) Then
Label26.Visible = True
p4 = "Yes"
Sarr(5) = 1
Else
Label26.Visible = False
p4 = "No"
End If
```

```
MSComm1.Output = "[2]"
Sleep (20)
Rtext = MSComm1.Input
pb = Val("&h" & Mid$(Rtext, 6, 2))
```

```
If ((pb And &H10) = &H10) Then
Label8.Visible = False
Saline.Visible = True
Sal = "Empty"
Sarr(6) = 1
```

```
Set objOutlookMsg = objOutlook.CreateItem(olMailItem)
With objOutlookMsg
.To = "9884169070@sms.sancharnet.in" ' The To Address Field
.Subject = "saline Emp" ' txtSubject - The Subject Field
.Body = "Saline bottle Empty"
.Importance = olImportanceHigh ' Priority. Mostly Not Required
.Send ' Send the Composed Message...
End With
```

```
Set objOutlookMsg = Nothing
End If
```

```
Else
Saline.Visible = False
Label8.Visible = True
Sal = "Not Empty"
End If
```

```
If (T1 > 100) Or (T2 > 100) Then
Label22.Visible = True
Else
Label22.Visible = False
End If
```

```

MSComm1.Output = "[2]"
Sleep (50)
Rtext = MSComm1.Input
pc = Val("&h" & Mid$(Rtext, 6, 2))
pc = pc And &HF

Select Case (pc And &HF)
Case &HE
    Atttxt(0) = "Rahul"
    Atttxt(1) = "Doctor"
    Atttxt(2) = Date
    Image1.Picture = LoadPicture("D:\PanimalarCBSS\dr1.gif")
    If vvv2 = 0 Then
        Atttxt(3) = Time
        Drs.AddNew
        Drs.Fields(0) = Atttxt(0)
        Drs.Fields(1) = Atttxt(1)
        Drs.Fields(2) = Atttxt(2)
        Drs.Fields(3) = Atttxt(3)
        Drs.Update
        vvv2 = 1
    End If
Case &HD
    Atttxt(0) = "Anita"
    Atttxt(1) = "Nurse"
    Atttxt(2) = Date
    Image1.Picture = LoadPicture("D:\PanimalarCBSS\dr2.gif")
    If vvv2 = 0 Then
        Atttxt(3) = Time
        Drs.AddNew
        Drs.Fields(0) = Atttxt(0)
        Drs.Fields(1) = Atttxt(1)
        Drs.Fields(2) = Atttxt(2)
        Drs.Fields(3) = Atttxt(3)
        Drs.Update
        vvv2 = 1
    End If
Case &HC
    Atttxt(0) = "Naresh"
    Atttxt(1) = "Compounder"
    Atttxt(2) = Date
    Image1.Picture = LoadPicture("D:\PanimalarCBSS\dr3.gif")
    If vvv2 = 0 Then
        Atttxt(3) = Time
        Drs.AddNew
        Drs.Fields(0) = Atttxt(0)

```



```

Drs.Fields(1) = Attxt(1)
Drs.Fields(2) = Attxt(2)
Drs.Fields(3) = Attxt(3)
Drs.Update
    vvv2 = 1
End If
Case &HB
    Attxt(0) = "Bhuvan"
    Attxt(1) = "Doctor"
    Attxt(2) = Date
    Image1.Picture = LoadPicture("D:\PanimalarCBSS\dr4.gif")
    If vvv2 = 0 Then
        Attxt(3) = Time
        Drs.AddNew
        Drs.Fields(0) = Attxt(0)
        Drs.Fields(1) = Attxt(1)
        Drs.Fields(2) = Attxt(2)
        Drs.Fields(3) = Attxt(3)
        Drs.Update
        vvv2 = 1
    End If
Case &HA
    Attxt(0) = "Raghu"
    Attxt(1) = "Doctor"
    Attxt(2) = Date
    Image1.Picture = LoadPicture("D:\PanimalarCBSS\dr1.gif")
    If vvv2 = 0 Then
        Attxt(3) = Time
        Drs.AddNew
        Drs.Fields(0) = Attxt(0)
        Drs.Fields(1) = Attxt(1)
        Drs.Fields(2) = Attxt(2)
        Drs.Fields(3) = Attxt(3)
        Drs.Update
        vvv2 = 1
    End If
Case &H9
    Attxt(0) = "Kamala"
    Attxt(1) = "Doctor"
    Attxt(2) = Date
    Image1.Picture = LoadPicture("D:\PanimalarCBSS\dr2.gif")
    If vvv2 = 0 Then
        Attxt(3) = Time
        Drs.AddNew
        Drs.Fields(0) = Attxt(0)
        Drs.Fields(1) = Attxt(1)
        Drs.Fields(2) = Attxt(2)

```

```

Drs.Fields(3) = Attxt(3)
Drs.Update
vvv2 = 1
End If
Case &H8
    Attxt(0) = "Soman"
    Attxt(1) = "Doctor"
    Attxt(2) = Date
    Image1.Picture = LoadPicture("D:\PanimalarCBSS\dr3.gif")
    If vvv2 = 0 Then
        Attxt(3) = Time
        Drs.AddNew
        Drs.Fields(0) = Attxt(0)
        Drs.Fields(1) = Attxt(1)
        Drs.Fields(2) = Attxt(2)
        Drs.Fields(3) = Attxt(3)
        Drs.Update
    vvv2 = 1
    End If
    Case &H7
        Attxt(0) = "Binu"
        Attxt(1) = "Doctor"
        Attxt(2) = Date
        Image1.Picture = LoadPicture("D:\PanimalarCBSS\dr4.gif")
        If vvv2 = 0 Then
            Attxt(3) = Time
            Drs.AddNew
            Drs.Fields(0) = Attxt(0)
            Drs.Fields(1) = Attxt(1)
            Drs.Fields(2) = Attxt(2)
            Drs.Fields(3) = Attxt(3)
            Drs.Update
            vvv2 = 1
        End If
    Case &H6
        Attxt(0) = "Algates"
        Attxt(1) = "Doctor"
        Attxt(2) = Date
        Image1.Picture = LoadPicture("D:\PanimalarCBSS\dr1.gif")
        If vvv2 = 0 Then
            Attxt(3) = Time
            Drs.AddNew
            Drs.Fields(0) = Attxt(0)
            Drs.Fields(1) = Attxt(1)
            Drs.Fields(2) = Attxt(2)
            Drs.Fields(3) = Attxt(3)
            Drs.Update

```

```

        vvv2 = 1
    End If
Case &H5
    Atttxt(0) = "Jagadeesh"
    Atttxt(1) = "Doctor"
    Atttxt(2) = Date
    Image1.Picture = LoadPicture("D:\PanimalarCBSS\dr2.gif")
    If vvv2 = 0 Then
        Atttxt(3) = Time
        Drs.AddNew
        Drs.Fields(0) = Atttxt(0)
        Drs.Fields(1) = Atttxt(1)
        Drs.Fields(2) = Atttxt(2)
        Drs.Fields(3) = Atttxt(3)
        Drs.Update
        vvv2 = 1
    End If
Case &H4
    Atttxt(0) = "Mohan"
    Atttxt(1) = "Wardboy"
    Atttxt(2) = Date
    Image1.Picture = LoadPicture("D:\PanimalarCBSS\dr3.gif")
    If vvv2 = 0 Then
        Atttxt(3) = Time
        Drs.AddNew
        Drs.Fields(0) = Atttxt(0)
        Drs.Fields(1) = Atttxt(1)
        Drs.Fields(2) = Atttxt(2)
        Drs.Fields(3) = Atttxt(3)
        Drs.Update
        vvv2 = 1
    End If
Case &H3
    Atttxt(0) = "Naga"
    Atttxt(1) = "Anaesthesian"
    Atttxt(2) = Date
    Image1.Picture = LoadPicture("D:\PanimalarCBSS\dr4.gif")
    If vvv2 = 0 Then
        Atttxt(3) = Time
        Drs.AddNew
        Drs.Fields(0) = Atttxt(0)
        Drs.Fields(1) = Atttxt(1)
        Drs.Fields(2) = Atttxt(2)
        Drs.Fields(3) = Atttxt(3)
        Drs.Update
        vvv2 = 1
    End If

```

Case &H2

Atttxt(0) = "Desai"

Atttxt(1) = "Nurse"

Atttxt(2) = Date

Image1.Picture = LoadPicture("D:\PanimalarCBSS\dr1.gif")

If vvv2 = 0 Then

Atttxt(3) = Time

Drs.AddNew

Drs.Fields(0) = Atttxt(0)

Drs.Fields(1) = Atttxt(1)

Drs.Fields(2) = Atttxt(2)

Drs.Fields(3) = Atttxt(3)

Drs.Update

vvv2 = 1

End If

Case &H1

Atttxt(0) = "Mukharjee"

Atttxt(1) = "Doctor"

Atttxt(2) = Date

Image1.Picture = LoadPicture("D:\PanimalarCBSS\dr2.gif")

If vvv2 = 0 Then

Atttxt(3) = Time

vvv2 = 1

End If

Case &HF

Atttxt(0) = ""

Atttxt(1) = ""

Atttxt(2) = ""

vvv2 = 0

Atttxt(3) = ""

Image1.Picture = LoadPicture("")

End Select

X = X + 1

If X >= 30 Then

rs.AddNew

rs(0) = Format(Date, "dd/mm/yyyy")

rs(1) = Format(Time, "hh:mm:ss")

rs(2) = T1 rs(3) = T2 rs(4) = Ecg1

rs.Update

Srs.AddNew

Srs(0) = Format(Date, "dd/mm/yyyy")

Srs(1) = Time Srs(2) = p1 Srs(3) = p2 Srs(4) = p3 Srs(5) = p4S rs(6) =Sal

Srs.Update

X = 0

End If

End Sub

Private Sub Timer2\_Timer()

Fcnt = Fcnt + 1

If Fcnt > 6 Then

Fcnt = 1

End If

If (Sarr(Fcnt) = 1) Then

MMControl1.FileName = Fnames(Fcnt)

MMControl1.Command = "Open"

MMControl1.Command = "Play"

Sarr(Fcnt) = 0

Timer2.Enabled = False

End If

End Sub

Private Sub Timer3\_Timer()

If Ecg1Timer.Enabled = True Then

If Strt = 1 And ptr > 0 Then

Egra1(1) = 50 - ((a(ptr) \* 10) / (200# / 50#))

ptr = ptr - 1

Else

Egra1(1) = 50 - (Ecg1 / (200# / 50#))

ptr = 19

Strt = 0

End If

For i = 100 To 2 Step -1

Egra1(i) = Egra1(i - 1)

Next

Picture1.Cls

vvv = 1

For i = 1 To 99

Picture1.Line (vvv, Egra1(i))-(vvv + 2, Egra1(i + 1)), vbYellow

vvv = vvv + 2

Next

EndIf

End Sub

Private Sub Timer4\_Timer()

```
If Heartbeattimer.Enabled = True Then
hcnt = hcnt + 1
Text2.Text = hcnt

If hcnt >= 60 Then
Text3.Text = Bcnt

If (Bcnt < 72 Or Bcnt > 84) Then
Label9.Caption = "Heart Beat Abnormal"
btimer.Enabled = True
Else
Label9.Caption = "Heart Beat Normal"
End If

hcnt = 0
Bcnt = 0
End If

If hcnt < 3 Then
Label9.Visible = True
Else
Label9.Visible = False
End If

End If

End Sub
```

**Form3.frm**

```
Private Sub Command1_Click()
```

```
Form3.Hide
```

```
Form5.Show
```

```
Form5.Data1.RecordSource="Doc"
```

```
Form5.Data1.Refresh
```

```
Form5.DBGrid1.ClearFields
```

```
Form5.DBGrid1.ReBind
```

```
End Sub
```

```
Private Sub Command2_Click()
```

```
Form3.Hide
```

```
Form5.Show
```

```
Form5.Data1.RecordSource="patient"
```

```
Form5.Data1.Refresh
```

```
Form5.DBGrid1.ClearFields
```

```
Form5.DBGrid1.ReBind
```

```
End Sub
```

```
Private Sub Command3_Click()
```

```
Form3.Hide
```

```
Form5.Show
```

```
Form5.Data1.RecordSource="Digital"
```

Form5.Data1.Refresh

Form5.DBGrid1.ClearFields

Form5.DBGrid1.ReBind

End Sub

Private Sub Command5\_Click()

Form3.Hide

Form5.Show

End Sub

Private Sub Command6\_Click()

End

End Sub

Private Sub Form\_DblClick()

Me.PrintForm

End Sub



Form5.frm

```
Private Sub Command1_Click()  
Data1.Recordset.MoveFirst  
While Not (Data1.Recordset.EOF)  
Data1.Recordset.Delete  
Data1.Recordset.MoveNext  
Wend  
DBGrid1.ClearFields  
End Sub
```

```
Private Sub Command2_Click()  
Form5.Hide  
Form4.Show  
Line(528,182)-(744,384),Form4.BackColor,BF  
oldx=536  
oldy=384  
k=0  
End Sub
```

```
Private Sub Command3_Click ()  
Form4.MSComm1.Output = "[3fff80ff]"  
Sleep(50)  
End  
End Sub
```

## Module1

### Option Explicit

Public db as Database

Public textval , dtext, SString, Sendstr, Rtext As String

Public trs, rps, prs, Srs, Drs As Recordset

Public X, Eout1, bout1, out, Timeout1, As Integer

Public pch, pa, pb, attin, pc, As Integer

Public oldx, oldy, Maxv As Single

Public k, grval As Integer

Public t1, tt0, tt1, t2, t3, t4, tt2, Ecg1 As Double

Public p1, p2, p3, p4, Sal As String

Public Strt As Integer

Public Fnames(1 to 7) As String

Public Egra1(1 to 100) As Single

Public Fcnt As Integer

Public a(0 to 19) As Single

Public pt1, pt2 As Single

Public ptr, vvv, vvv2 Sarr(1 to 7) As Integer

Public Beat, Bcnt, Bflg, Tcount As Single

Public hcnt, hout As Integer

Public Declare Function Sleep Lib "kernel32" (ByVal dwMilliseconds As Long) As Long

Public objOutlook As New Outlook.Application

Public objOutlookMsg As Outlook.MailItem

Public pflg, sale, t1flg, t2flg, hflg1, finflg, myflg As Integer

Public Lg As Long

# CONCLUSION

## **CONCLUSION**

This Project which demonstrates an automated patient monitoring system has its own merits and demerits which are discussed below:

### **Merits:**

2. The wireless alert system using WAP notifies physicians of critical results on their cellular phones.
3. With online recoding of medical parameters, the workload of the case providers and the nursing staff is reduced.
4. The clinical information database contains all data regarding the patients in electronic form.
5. The patient call switches help emergency situations to be handled quickly.

Future enhancements can be easily implemented with the PLC controller.

### **De-Merits:**

1. The heart beat sensor is highly temperature dependent and the dynamic characteristics change with different levels of ambient light and temperature level.
2. The ECG amplifier needs a high CMRR and additional narrowband filters are necessary to attenuate effects of mains and noise interference.
3. Thermistors used for respiratory temperature measurement may need to be heated initially in order to provide better sensitivity. However, this excessive heating may cause discomfort to the patient.

# **FUTURE ENHANCEMENTS**

## **FUTURE ENHANCEMENTS**

The human body scanning system could be made more sophisticated by incorporating blood pressure and EEG sensors. The analog channel inputs AN4 and AN7 can be used and the Port B can be programmed as an input port along with an additional ADC chip in the external circuit.

Hospitable –wide wireless capability would allows doctor to occur the patients' database using their word held computers.

The entire medical data acquisition could be made wireless and wearable. Such a package would contain the circuiting for inputs from ECG sensors, EEG sensors, pressure measurement and pulse rate transducers. This wearable module can transmit the data continuously over a fiber optic link or through an internet digital radio. The received data can be stored in separate memory and be processed by a microcontroller. This enhancement will enable monitoring of patients to be more flexible and strain-free.

# **COST ESTIMATE**

**COST ESTIMATE**

<b>COMPONENTS</b>	<b>PRICE (Rs.)</b>
Power supply board	105.00
Resistors and wires	35.00
Transistors	25.00
Trimpots	40.00
Switches	40.00
IC OP07	110.00
9 pin female connector	15.00
MAX232	70.00
PIC16F877	700.00
Heartbeat sensor	1300.00
ECG sensor	950.00
Saline bottle	17.00
PCB and GPB	50.00
Microcontroller board	800.00
Cabinet	450.00
Miscellaneous	600.00
<b>TOTAL</b>	<b>5307.00</b>



# **BIBLIOGRAPHY**

## **BIBLIOGRAPHY**

1.     **Handbook of Bio-Medical Instrumentation**     -     **R.S.Khandpur**
  
2.     **Bio-Medical Instrumentation**  
        **and measurements**                                 -     **Leslie Cromwell**  
    -     **Fred. J. Wejnbell**  
    -     **Erich . A. Pleiffer**
  
3.     **Linear Integrated Circuits**                                 -     **Roy Chowdary**
  
4.     **IBM PC Handbook**   -     **IBM Corporation**
  
5.     **www.microchip.com**
  
6.     **VISUAL BASIC---- Ground Up**                                 -     **Gary Cornell**